



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-006691</b>	Instrument Location <b>KENMARE PD</b>
Reason for Install/Repair	
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <b>Annual Inspection at Crime Laboratory</b>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): **N/A**
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <b>P. Hall #133223</b>	Date <b>5-11-23</b>
---	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <b>Ame Nahl</b>	Certified Date <b>12May2023</b>
--	------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006691  
Location = KMPD                    8164.14.00 09/16  
05/11/2023                            16:33

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\* () \_ + ?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\* () \_ + ?

Current Instrument Setup

Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer (Y/N)?                No  
Display Volume?                      No  
Disable On Memfull?                Yes  
# of Print Copies?                   1  
Select Std (D/W/I)?                Dry  
Standard Value?                      0.080  
Standard Lot #?                      01923080A3  
Standard Cyl #?                      61  
Standard Expiration?                02/05/2025  
Oper No?                              133223

Flow Cal. Date:                      08/28/2018  
Slope                                  699  
Intercept                              -660610

IR Calibration Date:                08/28/2018  
    3um                                  9um

-----  
0th Coef (\*100):                    -9426                                  -20764  
1st Coef (\*100):                    258193                                134855  
2nd Coef (\*100):                    3532                                    1465  
H2O adj (mg/l\*10k):                408                                      547

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ALLEN HALL

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006691  
Location = KMPD                    8164.14.00 09/16  
05/11/2023                            16:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:35
02 Std. Gas	0.080	16:35
03 Room Air	0.000	16:36
04 Std. Gas	0.080	16:36
05 Room Air	0.000	16:37
06 Std. Gas	0.080	16:37
07 Room Air	0.000	16:38

Lot No = 01923080A3  
Cyl No = 61  
Exp Date = 02/05/2025  
County = 51                            Oper No. = 133223



133223

Operator Signature  
ALLEN HALL

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006691  
Location = KMPD      8164.14.00 09/16  
05/11/2023      16:39

Test	AC	Time
01 Room Air	0.000	16:41
02 *Subject Test	RFI*	16:42
03 Room Air	0.000	16:42

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 05/11/2000  
Sub Sex = Unknown      Weight = NA  
Test = OTH      Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 01923080A3  
Cyl No = 61  
Expiration Date = 02/05/2025  
County = 51      Oper No. = 133223

  
Operator Signature  
ALLEN HALL

Remarks:

Form 106-I8000