



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 5937

A. General Setup and Checks:

1. Initial start up
  - Power On
  - 5-minute countdown during initial warm-up
  - 1-minute countdown to begin test after standby mode
2.  Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
  - Level 1 (Operator)
  - Level 2 (Field Inspector)
  - Level 3 (Lab Personnel)
4.  Printer Setup (Level 3,S,P)
5.  Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1.  Internal Printer Test (Level 1,P). Sign and attach record.
  - Paper Feed
2.  External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
  - Operator ID
  - Invalid ID
  - ND License
  - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1.  Continuous Air Blank (Level 2,A)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 1 of 4

2.  Diagnostic (Level 2,D). Sign and attach test record.

3.  Prelim Data Entry (Level 2,E)

4. Ser Number/Version Info (Level 2,V):

Inst Ser#: 80-005937 DSP Ver: 1362.14  
Core Ver: 816f.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 13170 s: 2 C: 46.93  
9: 13528 s: 1 B: 44

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0→1

2. Function (F) Keys:

F5 – Tone                       F9 - Fan  
 F6 – Green LED                 F10 – Pump  
 F7 – Red LED                     F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 954 CF: 1.06

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 203 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 44.5 I: 0.0  
B: 45.0 V: 0.0

Cell Temperature (C) between  $47 \pm 0.02$  °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 468 psi Regulator: 480 psi

Display and Regulator  $\pm 50$  psi of each other  Yes or  No

Gas tank tare necessary?  Yes or  No

If Yes, display readings after tare (Level 3,M,C,G):

Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

7.  General Setup (Level 3,S,G)

8.  Configure Start Test (Level 3,S,T)

9.  Network (Level 3,S,N)

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 2 of 4

E. Tests (Sign and attach all test records):

1.  Configure wet bath simulator for the following test (Level 1,S)  
 Zero AC-Wet Calibration Check (Level 1,C):  
 Sim ser #: MP3066 Lot #: 33023 Exp date: N/A  
 Known value: 0.000 AC  
 Results are Zero AC
2.  Configure wet bath simulator for the following test (Level 1,S)  
 Low AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\leq 0.03$  AC: 0.020 AC  
 Sim ser #: MPS318 Lot #: 202201A Exp date: 11Jan24  
 Results  $\pm 0.005$  of known AC
3.  Configure wet bath simulator for the following test (Level 1,S)  
 High AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\geq 0.25$  AC: 0.300 AC  
 Sim ser #: MP3067 Lot #: 202201F Exp date: 18Jan24  
 Results  $\pm 5\%$  of known AC
4.  Configure dry gas standard for the following tests (Level 1,S)  
 Known Value: 0.080 AC  
 Gas Cylinder Lot #: 26021080A1 Cylinder #: 25 Exp. Date: 10/5/23
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7345 Lot #: 1057 Exp. Date: N/A  
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7349 Lot #: ACE3 Exp. Date: N/A  
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):  
 Sim. Ser #: DR7347 Lot #: MET3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):  
 Sim. Ser #: DR7350 Lot #: IPA3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
9.  Ambient Fail (Level 1,B): Display read and test record correct.
10.  Invalid Sample (Level 1,B): Display read and test record correct.
11.  Sequence Aborted (Level 1,B): Display read test record correct.
12.  ABA (Level 1,B): Test record correct.

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 3 of 4

- 13.  Calibration check out of tolerance (CMS mode): Display read and test record correct.
- 14.  RFI Detected (B,C, or CMS mode): Display read and test record correct.
- 15.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
- 16.  Deficient Samples (CMS mode): Display read and test record correct.
- 17.  Improper Sample (CMS mode): Display read and test record correct.
- 18.  Difference Too Great (CMS mode): Display read and test record correct.
- 19.  CMS (CMS mode): Test record correct.
- 20. Stability Test (Level 3,T,S):

Ave: 0.0776 Std Dev: 0.005 Rel Std Dev: 0.65%  
*0.0005  
15 May 23  
Jaw*

Maintenance/Repairs/Notes: Replaced modem port. Had to open  
Intox 50 then performed Initial Inspection.

Instrument is acceptable to be used in the field.  Yes or  No

If No, state reason(s) why:

*Am Mel*  
 Inspector Signature

12 May 2023  
 Date

*Janelle Pertschiller*  
 Reviewer

15 May 2023  
 Date

Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11784 Revision: 2
Intoxilyzer 8000 Initial Inspection	Status: Published Published: 04/24/2023
UNCONTROLLED WHEN PRINTED	Page 4 of 4



Intoxilyzer Test Record and Checklist  
NDDAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-005937  
Location = TOXL 8164.14.00 09/16  
05/12/2023 09:43

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Serial: 34921080A1  
Data Entry Mode: Standard  
Start Test Sequence: ABCACABA  
Display Prelim Result: Yes  
Display Third Result: Yes  
Inhib Printer: No  
Display Volume: No  
Disable On MemFull?: Yes  
# of Print Copies: 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 34921080A1  
Standard Cyl #? 55  
Standard Expiration? 02/05/2024  
Oper No? 133237

Flow Cal. Date: 11/17/2017  
Slope: 606  
Intercept: -534563

IR Calibration Date: 11/17/2017  
3um 9um


0th Coef(\*100): -30534 -23772

1st Coef(\*100): 294336 142357

2nd Coef(\*100): -906 850

H2O adj(ng/l\*10K): 205 250

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks:

Internal Print Test

Form 106-18000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                8164.14.00 09/16  
05/12/2023                        09:44

Test	AC	Time
01 Room Air	0.000	09:45
02 *Subject Test	REF*	09:45
03 Room Air	0.000	09:46

\*Subject Test Refused

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = REFUSAL  
Dr. Lic. = ND/TES989643  
Lot No = 34921080A1  
Cyl No = 55  
Expiration Date = 02/05/2024  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Refusal

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                 8164.14.00 09/16  
05/12/2023                        09:43

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks:

Diagnostic

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005937  
Location = TOXL      8164.14.00 09/16  
05/12/2023      09:51

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:51
02 Std. Sol.	0.000	09:52
03 Room Air	0.000	09:53
04 Std. Sol.	0.000	09:53
05 Room Air	0.000	09:54
06 Std. Sol.	0.000	09:54
07 Room Air	0.000	09:55

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 33023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Zero Ac Check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005937  
Location = TOXL      8164.14.00 09/16  
05/12/2023      09:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:57
02 Std. Sol.	0.021	09:58
03 Room Air	0.000	09:59
04 Std. Sol.	0.020	09:59
05 Room Air	0.000	10:00
06 Std. Sol.	0.021	10:01
07 Room Air	0.000	10:01

08 Sim Temp = 34.0°C

Simul Ser No = MP5318

Std Sol No = 202201A

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks:

LOW AC Check - 0.020AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005937  
Location = TOXL      8164.14.00 09/16  
05/12/2023      10:02

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:02
02 Std. Sol.	0.301	10:03
03 Room Air	0.000	10:04
04 Std. Sol.	0.302	10:04
05 Room Air	0.000	10:05
06 Std. Sol.	0.305	10:06
07 Room Air	0.000	10:06


08 Sim Temp = 34.0°C

Simul Ser No = MP3067

Std Sol No = 202201F

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: High Ac Check - 0.300Ac

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            10:12

Test	AC	Time
01 Room Air	0.000	10:12
02 *Subject Test	INT*	10:13
03 Room Air	0.000	10:14

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent Check  
              ETOH + Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            10:14

Test	AC	Time
01 Room Air	0.000	10:15
02 *Subject Test	INT*	10:15
03 Room Air	0.000	10:16

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                            Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Interferent check  
Acetone*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            10:16

Test	AC	Time
01 Room Air	0.000	10:17
02 *Subject Test	INT*	10:18
03 Room Air	0.000	10:18

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
Methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005937  
Location = TOXL      8164.14.00 09/16  
05/12/2023      10:19

Test	AC	Time
01 Room Air	0.000	10:19
02 *Subject Test	INT*	10:20
03 Room Air	0.000	10:20

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male      Weight = NA  
Test = OTH      Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check  
IPA

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                8164.14.00 09/16  
05/12/2023                        10:20

Test	AC	Time
01 Room Air	AMB*	10:21
02 Room Air	0.000	10:22

\*Invalid Test  
Check Ambient Conditions

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = AMBIENT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *IS Ambient Fail*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                8164.14.00 09/16  
05/12/2023                        10:24

Test	AC	Time
01 Room Air	0.000	10:25
02 Invalid Sample	X.XXX	10:25
03 Room Air	0.000	10:26

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = INVALID  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Invalid Sample*

Form 106-I8000






Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                 8164.14.00 09/16  
05/12/2023                        10:32

Test	AC	Time
01 Room Air	0.000	10:32
02 Subject Test 1	0.000	10:33
03 Room Air	0.000	10:34
04 Reported AC	0.000	10:33

No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                     Weight = NA  
Test = OTH                         Cit = ABA  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: ABA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                8164.14.00 09/16  
05/12/2023                      10:35

Test	AC	Time
01 Diagnostic	OK	10:36
02 Room Air	0.000	10:37
03 Subject Test 1	0.000	10:37
04 Room Air	0.000	10:39
05 Std. Gas	0.047*	10:40
06 Room Air	0.000	10:41

\*Cal Check Out of Tolerance

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = CCOOT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                          Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Cal check out of tolerance

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                8164.14.00 09/16  
05/12/2023                        10:44

Test	AC	Time
01 Room Air	0.000	10:45
02 *Subject Test	RFI*	10:45
03 Room Air	0.000	10:45

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = RFI CHECK  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *RFI Check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            10:47

Test	AC	Time
01 Diagnostic	OK	10:48
02 Room Air	0.000	10:49
03 *Subject Test	0.000*	10:52
04 Room Air	0.000	10:54
05 Std. Gas	0.077	10:55
06 Room Air	0.000	10:56
07 *Subject Test	0.000*	11:00
08 Room Air	0.000	11:01

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = DEFICIENT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Deficient Sample*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005937  
Location = TOXL                 8164.14.00 09/16  
05/12/2023                        11:04

Test	AC	Time
01 Diagnostic	OK	11:05
02 Room Air	0.000	11:06
03 Subject Test 1	0.000	11:06
04 Room Air	0.000	11:08
05 Reference	IPS*	11:08
06 Room Air	0.000	11:09

\*Improper Sample

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                     Weight = NA  
Test = OTH                         Cit = IMPROPER  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                        Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Improper Sample*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            11:09

Test	AC	Time
01 Diagnostic	OK	11:10
02 Room Air	0.000	11:10
03 Subject Test 1	0.041*	11:11
04 Room Air	0.000	11:13
05 Std. Gas	0.078	11:14
06 Room Air	0.000	11:15
07 Subject Test 2	0.000*	11:16
08 Room Air	0.000	11:17

\*Difference Too Great

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                            Cit = DIFF TOO GREAT  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                            Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005937  
Location = TOXL                    8164.14.00 09/16  
05/12/2023                            11:18

Test	AC	Time
01 Diagnostic	OK	11:18
02 Room Air	0.000	11:19
03 Subject Test 1	0.000	11:19
04 Room Air	0.000	11:22
05 Std. Gas	0.077	11:23
06 Room Air	0.000	11:25
07 Subject Test 2	0.000	11:25
08 Room Air	0.000	11:26
09 Reported AC	0.000	11:19

Difference OK  
No RFI Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                            Weight = NA  
Test = OTH                                Cit = CMS TEST  
Dr. Lic. = ND/TES989643  
Lot No = 26021080A1  
Cyl No = 25  
Expiration Date = 10/05/2023  
County = 08                                Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: CMS Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005937  
Location = TOXL      8164.14.00 09/16  
05/12/2023      11:27

Test	AC	Time
01 Room Air	0.000	11:28
02 Std. Gas	0.078	11:28
03 Room Air	0.000	11:29
04 Std. Gas	0.077	11:29
05 Room Air	0.000	11:30
06 Std. Gas	0.078	11:30
07 Room Air	0.000	11:31
08 Std. Gas	0.078	11:31
09 Room Air	0.000	11:32
10 Std. Gas	0.078	11:32
11 Room Air	0.000	11:32
12 Std. Gas	0.078	11:33
13 Room Air	0.000	11:33
14 Std. Gas	0.078	11:34
15 Room Air	0.000	11:34
16 Std. Gas	0.078	11:35
17 Room Air	0.000	11:35
18 Std. Gas	0.077	11:36
19 Room Air	0.000	11:36
20 Std. Gas	0.077	11:36
21 Room Air	0.000	11:37
22 Std. Gas	0.078	11:37
23 Room Air	0.000	11:38
24 Std. Gas	0.078	11:38
25 Room Air	0.000	11:39
26 Std. Gas	0.077	11:39
27 Room Air	0.000	11:40
28 Std. Gas	0.077	11:40
29 Room Air	0.000	11:41
30 Std. Gas	0.077	11:41
31 Room Air	0.000	11:42
32 Std. Gas	0.077	11:42
33 Room Air	0.000	11:42
34 Std. Gas	0.078	11:43
35 Room Air	0.000	11:43
36 Std. Gas	0.077	11:44
37 Room Air	0.000	11:44
38 Std. Gas	0.078	11:45
39 Room Air	0.000	11:45
40 Std. Gas	0.078	11:46
41 Room Air	0.000	11:46
42 Cal Check Stats:		
43 Average	0.0776	
44 Std Dev	0.0005	
45 Rel Std Dev	0.65%	



Operator Signature  
ANNA NAREHOOD

Remarks:

*On hold*  
*Stability Test*  
Form 106-18000