

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

Serial Number 80 - 005952	Instrument Location Wisher P.	D			
Reason for Install/Repair					
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change				
Other (Specify) Annual Inspection at Co	rime Laboratory				
	<i>d</i>				
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or	C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Age	ncy.			
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form			
Field Inspector Signature		Date 5 13-23			
Crime Laboratory Use Only This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on	file at the Office of Attorney General, Offi	rmine alcohol concentration from ne Laboratory Division, in the			
County of Burleigh, North Dakota, is certified to be a true and correct c	opy of the documents received.				

Reviewed/Certified By

Certified Date

15May2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005952 Location = WISH 8164.14.00 09/16 05/13/2023 22:18

******* Printer Test ********

abcdefqhijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABA Display Prelim Rslt? Yes DABACABA

Display Third Digit? Yes Inhib Printer(Y/N)? No No Yes Display Volume? Disable On Memfull? # of Print Copies? 1 # OI Print Copies.
Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 34921080A1

38 Standard Expiration? 02/05/2024

Oper No? 010912

Flow Cal. Date: 04/16/2021 Slope 657

IR Calibration Date: 04/16/2021

3um 9um

-744615

Oth Coef(*100): -24450 1st Coef(*100): 266154 2nd Coef(*100): 2284 -14639 132500 1091 H20 adj(mq/l*10k): 516323

****** Printer Test End *******

Operator Signature ANDREW TUNHEIM

Remarks:

Intercept

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WISH
05/13/2023

Alcohol Analyzer SN 80-005952 8164.14.00 09/16 22:19

DRY CAL CHECK

Test			AC	Time	
	01	Room	Air	0.000	22:20
	02	Std.	Gas	0.080	22:20
	03	Room	Air	0.000	22:21
	04	Std.	Gas	0.080	22:21
	05	Room	Air	0.000	22:22
	06	Std.	Gas	0.080	22:22
	07	Room	Air	0.000	22:23

Lot No = 34921080A1

Cyl No = 38

Exp Date = 02/05/2024

County = 26

Oper No. = 010912

Operator Signature ANDREW TUNHEIM

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WISH
05/13/2023

Alcohol Analyzer SN 80-005952 8164.14.00 09/16 22:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:29
02 Std. Gas	0.080	22:29
03 Room Air	0.000	22:29
04 Reference	RFI*	22:30
05 Room Air	0.000	22:30

*Invalid Test Inhibited - RFI

Lot No = 34921080A1

Cyl No = 38

Exp Date = 02/05/2024

County = 26

Oper No. = 010912

Operator Signature ANDREW TUNHEIM

Remarks:

Form 106-I8000