



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59281 (06/2018)

| | |
|---|------------------------------------|
| Serial Number 80-007091 | Instrument Location Mott |
| Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Annual Inspection at crime Laboratory</u> | |

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

| | |
|---|--------------------------|
| Field Inspector Signature <i>Sarah D. Abramo</i> | Date 5/15/2023 |
|---|--------------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | |
|--|------------------------------------|
| Reviewed/Certified By <i>Ame Nahl</i> | Certified Date 16May2023 |
|--|------------------------------------|

Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

| | |
|-------------------------|------------------|
| CMI, Inc. Intoxilyzer | Alcohol Analyzer |
| North Dakota Model 8000 | SN 80-007091 |
| Location = MOTT | 8164.16.00 09/18 |
| 05/15/2023 | 10:36 |

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^&*()_+?

Current Instrument Setup


| | |
|----------------------|------------|
| Data Entry Mode: | Enabled |
| Start Test Sequence: | DABACABA |
| Display Prelim Rslt? | Yes |
| Display Third Digit? | Yes |
| Inhib Printer(Y/N)? | No |
| Display Volume? | No |
| Disable On Memfull? | Yes |
| # of Print Copies? | 1 |
| Select Std (D/W/I)? | Dry |
| Standard Value? | 0.080 |
| Standard Lot #? | 26021080A1 |
| Standard Cyl #? | 9 |
| Standard Expiration? | 10/05/2023 |
| Oper No? | 100620 |

| | |
|-----------------|------------|
| Flow Cal. Date: | 09/25/2018 |
| Slope | 681 |
| Intercept | -651030 |

| | |
|----------------------|------------|
| IR Calibration Date: | 09/25/2018 |
| 3um | 9um |

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -39121 | -20906 |
| 1st Coef(*100): | 240667 | 132345 |
| 2nd Coef(*100): | 3960 | 1516 |
| H2O adj(mg/l*10k): | 899 | 428 |

***** Printer Test End *****



Operator Signature
 SARAH WARNER

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007091
Location = MOTT 8164.16.00 09/18
05/15/2023 10:41

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:42 |
| 02 Std. Gas | 0.078 | 10:42 |
| 03 Room Air | 0.000 | 10:43 |
| 04 Std. Gas | 0.078 | 10:43 |
| 05 Room Air | 0.000 | 10:43 |
| 06 Std. Gas | 0.077 | 10:44 |
| 07 Room Air | 0.000 | 10:44 |

Lot No = 26021080A1
Cyl No = 9
Exp Date = 10/05/2023
County = 21 Oper No. = 100620


Operator Signature
SARAH WARNER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007091
Location = MOTT 8164.16.00 09/18
05/15/2023 10:45

DRY CAL CHECK

| Test | AC | Time |
|-------------|------|-------|
| 01 Room Air | RFI* | 10:46 |
| 02 Room Air | RFI* | 10:46 |

*Invalid Test
Inhibited - RFI

Lot No = 26021080A1
Cyl No = 9
Exp Date = 10/05/2023
County = 21 Oper No. = 100620



Operator Signature
SARAH WARNER

Remarks:

Form 106-I8000