

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 4194	Instrument Location Burleigh Moton 1	Detention
Reason for Install/Repair		•
☐ Install After Receiving From Crime Laboratory ☐ Install Aft	er Location Change	
Dother (Specify) Annual Inspection a	+ Crime Laboratory	
Check When Done:	O	
<ol> <li>Surge Protector Installed/Property Grounded.</li> </ol>		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Local	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C). ACA test is a	Hached- AEN	
C. Radio Frequency Interference (RFI) Test (CMS Mode of		ring Test)
7. Repair and/or Maintenance Performed (if any): Non-	e Nueded	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN:	50496, Form 120-G) and Place it by the In	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Location at the Age	ncy.
<ol> <li>Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.</li> </ol>	ntoxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date , ,
	809	5/15/23
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct of	file at the Office of Attorney General, Crin	
Reviewed/Certified By		Certified Date
		221/00112

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004194
Location = BMDC	8164.14.00 09/16
05/15/2023	10:24

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current In	strument	Setup
------------	----------	-------

earrene riberament becap	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	01923080A3
Standard Cyl #?	53
Standard Expiration?	02/05/2025
Oper No?	130315

Flow Cal.	Date:	06/23/2009
Slope		681

Intercept -622430

IR	Calibration	Date:	05/04/2016
		3um	9um

	· -			
0th	Coef(*100);		-17808	-24119
1st	Coef(*100):		263893	132021
2nd	Coef(*100):		3080	1296
H20	adj(mg/l*10k	):	403	519

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature JARED LEMIEUX

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
05/15/2023

Alcohol Analyzer SN 80-004194 8164.14.00 09/16 10:27

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	10:27
02	std.	Gas	0.081	10:27
03	Room	Air	0.000	10:28
04	Std.	Gas	0.080	10:28
05	Room	Air	0.000	10:29
06	std.	Gas	0.080	10:29
07	Room	Air	0.000	10:30

Lot No = 01923080A3

Cyl No = 53

Exp Date = 02/05/2025

County = 08

Oper No. = 130315

Operator Signature JARED LEMIEUX

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BMDC
05/15/2023

Alcohol Analyzer SN 80-004194 8164.14.00 09/16 10:25

## DRY CAL CHECK

Test			AC	Time	
01	Room	Air	0.000	10:25	
02	Std.	Gas	0.080	10:26	
03	Room	Air	RFI*	10:26	
04	Room	Air	RFI*	10:26	

\*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 53

Exp Date = 02/05/2025

County = 08

Oper No. = 130315

Operator Signature
JARED LEMIEUX

Remarks:

Form 106-I8000