

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 004300	Instrument Location NAPO		
Reason for Install/Repair			
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	er Location Change		
Other (Specify) Annual Inspection at	crime Laboratory		
Check When Done:	0		
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3 Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
∡ 6. Run Tests:			
A. Print Test (Level1, Function P).		,	
☑ B. ACA Test (Level 1, Function C).			
💢 C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).	
∑_7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.	
🥱 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.	
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature Michael Walter		Date 5 - 22 - 23	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fil County of Burleigh, North Dakota, is certified to be a true and correct co	e at the Office of Attorney General, Crin		
Reviewed/Certified By And Mol		Certified Date 23May2023	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = NAPO 8164.14.00 09/16
05/22/2023 14:15

****** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Dry Select Std (D/W/I)? 0.080 Standard Value? Standard Lot #? 34921080A1 Standard Cyl #? 59 Standard Expiration? 02/05/2024 131100 Oper No?

Flow Cal. Date: 06/25/2019 Slope 614 Intercept -551153

IR Calibration Date: 06/25/2019
3um 9um

 0th Coef(*100):
 -18792
 -22164

 1st Coef(*100):
 267001
 140777

 2nd Coef(*100):
 1903
 1083

 H2O adj(mg/l*10k):
 403
 445

****** Printer Test End *******

Operator Signature MICHAEL MAROLT

1111

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = NAPO 05/22/2023 SN 80-004200 8164.14.00 09/16 14:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:18
02 Std. Gas	0.080	14:18
03 Room Air	0.000	14:19
04 Std. Gas	0.080	14:19
05 Room Air	0.000	14:20
06 Std. Gas	0.080	14:20
07 Room Air	0.000	14:20

Lot No = 34921080A1

Cyl No = 59

Exp Date = 02/05/2024

County = 24

Oper No. = 131100

Operator Signature MICHAEL MAROLT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-004200 Location = NAPO 05/22/2023

Alcohol Analyzer 8164.14.00 09/16 14:25

Test		AC	Time
	Room Room	 RFI* 0.000	14:28 14:28

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 07/20/1999

Sub Sex = MaleWeight = 600Cit = 9999999Test = DUI

Dr. Lic. = ND/TES946666

Lot No = 34921080A1

Cyl No = 59

Expiration Date = 02/05/2024

County = 24

Oper No. = 131100

Operator Signature MICHAEL MAROLT

Remarks:

Form 106-I8000