	NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION								
	INTOXILYZER® 8000 ANNUAL II								
Into	Intoxilyzer [®] 8000 Serial Number: 80-00 5944 Inspection Location: ToxL								
Α.	 Pre-Inspection 1. Items with Instrument: Gas Cylinder Ves or No (If Yes, Lot #_ Keys Yes or O Power Cord Yes or No 2. Download Data 3. Upload Operator File 4. Current Location Code: <u>CASS</u> 5. Eattery Check Was the external battery pack replace 6. O-Rings Replaced Simulator O-Ring Replaced Breath Tube O-Ring 	20021080 A1 Cyl. # 20) d? Yes or No Yes or No Yes or No							
В.	 General Setup and Checks: 1. X Diagnostics passed and instrument in "Re 2. A Breath tube heated 3. X Date, time and location code (Level 2,E). Time Zone: CST or MDT (Time on test 4. X Print test (Level 1,P). Sign and attach test 5. X Tank monitor (Level 3,D,G). Display:psi Regulator: Display and Regulator ± 50 psi of each Gas tank tare necessary? If Yes, display readings after tare (Lev Display:psi Regulator: 	Re-set if necessary. st records will be in time zone circled). t record. <u>Solution</u> psi n other (resor No Yes or No el 3,M,C,G):							
C.	 Tests (Sign and attach test records): 1. ☆ Configure simulator for the following test of Wet Calibration Check - Low AC (Level Known Value ≤ 0.03 AC: 0.020 Sim. Ser #:	el 1,C) _AC 							
Тох	icology Unit - Breath Alcohol Section	Qualtrax ID: 11698 Revision: 1 Status: Published							
Into	xilyzer 8000 Annual Inspection	Published: 03/24/2023							
UN	CONTROLLED WHEN PRINTED	Page 1 of 2							

2. X Configure simulator for the follow Wet Calibration Check - High							
Known Value ≥ 0.25 AC: _0	300 AC						
Sim. Ser #:	e7						
Lot #: <u>20ス入り</u> Exp. Date: <u>18JCu</u>							
Results ± 5% AC of known							
3. Configure dry gas standard for th							
Known Value : 0.000 A	NC .						
Gas Cylinder Lot #:	21080A1						
Cylinder #:24							
Exp. Date: 10/5/23							
4. A Interferent Check (Level 1,B)	5% Acotopo						
Known Value: <u>0.10</u> AC + 0.05 Sim. Ser #: <u>DR7345</u>	5% Acelone						
Lot #: ICS 7							
Exp. Date: N(A							
Display reads "Interferent	Detect"						
5. 🔽 RFI Check (CMS Mode)							
Display reads "RFI Detect	"						
6. 💢 Dry Calibration Check (Level 1,C							
Test 1 0.081 Test 4 0							
Test 2 <u>0.081</u> Test 5 <u>(</u>	.081 Test 8 0.081						
Test 3 <u>0.080</u> Test 6 (
Average							
🔀 Results ± 0.005 AC of knc	wh AC						
D. Remarks/Maintenance: Added plech	ical tape around whip rip covering						
in mouth breathtube. The brai	ided covering was trayed.						
	5 1						
/	3						
Instrument is acceptable to be used in the field.	res or No						
If No. state reason(s) why:							
If No, state reason(s) why:							
If Yes, change location code back to A.4.							
	25 May 2023						
nspestor Signature Date							
Sanelle Pertschuller	25114122022						
Reviewer	Date						
	Duio						
Toxicology Unit - Breath Alcohol Section	Qualtrax ID: 11698 Revision: 1						
Texteelegy entry breath Alcohol Section	Status: Published						
Intoxilyzer 8000 Annual Inspection	Published: 03/24/2023						
Page 2 of							

UNCONTROLLED WHEN PRINTED

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005944
Location = TOXL	8164.14.00 09/16
05/25/2023	10:59

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	26021080A1
Standard Cyl #?	20
Standard Expiration?	10/05/2023
Oper No?	133237
Flow Cal. Date:	06/20/2013
Slope	677
Intercept	-702620

IR Calibration Date: 06/09/2017 3um 9um Oth Coef(*100): -16158 -17815 1st Coef(*100): 266780 136962 2nd Coef(*100): 1117 1089 H2O adj(mg/l*10k): 317 364

********** Printer Test End **********

Operator Signature ANNA NAREHOOD

Remarks:

not Test

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005944Location = TOXL8164.14.00 09/1605/25/202311:21

	WET CAL CHECK	Ĺ		
Time	AC		est	Τe
11:22	0.000	Air	Room	01
11:22	0.021	Sol.	Std.	02
11:23	0.000	Air	Room	03
11:24	0.021	Sol.	Std.	04
11:24	0.000	Air	Room	05
11:25	0.021	Sol.	Std.	06
11:25	0.000	Air	Room	07

08 Sim Temp = 34.0°C

Simul Ser No = MP5318 Std Sol No = 202201A County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

LOW ACCHELK-0.020AC

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005944Location = TOXL8164.14.00 09/1605/25/202311:31

			WET CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	11:32
02	Std.	Sol.	0.304	11:33
03	Room	Air	0.000	11:33
04	Std.	Sol.	0.305	11:34
05	Room	Air	0.000	11:35
06	Std.	Sol.	0.304	11:35
07	Room	Air	0.000	11:36

08 Sim Temp = 34.0°C

Operator Signature ANNA NAREHOOD

Remarks:

High Ac Check- 0.300Ac

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005944 Location = TOXL 8164.14.00 09/16 05/25/2023 11:37 AC Time Test 01 Room Air 0.000 11:38 02 *Subject Test INT* 11:38 03 Room Air 0.000 11:39 *Invalid Test Interferent Detected Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Test = OTHCit = INTERFERENT CK Dr. Lic. = ND/TES989643 Lot No = 26021080A1Cyl No = 24Expiration Date = 10/05/2023County = 08Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Interforent Cheur

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005944 8164.14.00 09/16 Location = TOXL11:39 05/25/2023 Time AC Test 11:40 OK 01 Diagnostic 11:40 02 Room Air 0.000 11:40 03 *Subject Test RFI* 11:41 04 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Weight = NA Sub Sex = Male Cit = RFI CHECK Test = OTHDr. Lic. = ND/TES989643 Lot No = 26021080A1Cyl No = 24Expiration Date = 10/05/2023Oper No. = 133237 County = 08

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks:

Cherr

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005944Location = TOXL8164.14.00 09/1605/25/202311:41

			DRY CAL CHECK	
Т	est		AC	Time
01	Room	Air	0.000	11:42
02	Std.	Gas	0.081	11:42
03	Room	Air	0.000	11:42
04	Std.	Gas	0.081	11:43
05	Room	Air	0.000	11:43
06	Std.	Gas	0.080	11:44
07	Room	Air	0.000	11:44

Lot No = 26021080A1 Cyl No = 24 Exp Date = 10/05/2023 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Dry calcherk #1-3

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005944Location = TOXL8164.14.00 09/1605/25/202311:44

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	11:45
02	Std.	Gas	0.081	11:45
03	Room	Air	0.000	11:46
04	Std.	Gas	0.081	11:46
05	Room	Air	0.000	11:47
06	Std.	Gas	0.082	11:47
07	Room	Air	0.000	11:48

Lot No = 26021080A1 Cyl No = 24 Exp Date = 10/05/2023 County = 08 Op

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Check #446 U

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005944Location = TOXL8164.14.00 09/1605/25/202311:48

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	11:49
02	Std.	Gas	0.081	11:49
03	Room	Air	0.000	11:49
04	Std.	Gas	0.081	11:50
05	Room	Air	0.000	11:50
06	Std.	Gas	0.082	11:51
07	Room	Air	0.000	11:51

Lot No = 26021080A1 Cyl No = 24 Exp Date = 10/05/2023 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Cal chere #7-9

Remarks: