

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006678 Instrument Location 6F	PD
Reason for Install/Repair	
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change	
Dother (Specify) Annual Inspection at Crime Laboratory	
THE STATE OF THE S	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
🛱 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Do	uring Test).
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair ( 104-G), Print Test, ACA Test, and RFI Test.	Checkout (SFN59281, Fðrm
Field Inspector Signature Lawrei	Date 05-26-23
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to dete the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Cri County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	
Reviewed/Certified By	Certified Date 26May 2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = GFPD 8164.14.00 09/16
05/26/2023 08:26

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? NoDisplay Volume? No Disable On Memfull? Yes # of Print Copies? 1 Dry Select Std (D/W/I)? Standard Value? 0.080 01923080A3 Standard Lot #? Standard Cyl #? 57 Standard Expiration? 02/05/2025 039947 Oper No?

Flow Cal. Date: 08/22/2016 Slope 686

Intercept -611596

IR Calibration Date: 08/22/2016
3um 9um

Oth Coef(\*100): -15445 -21650

1st Coef(\*100): 251854 131766

2nd Coef(\*100): 3094 1378

H2O adj(mg/l\*10k): 549 513

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
JONATHAN LAMPI

Remarks:

Annual Inspection Install
Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = GFPD 05/26/2023

SN 80-006678 8164.14.00 09/16 08:44

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:45
02 Std. Gas	0.081	08:45
03 Room Air	0.000	08:46
04 Std. Gas	0.081	08:46
05 Room Air	0.000	08:47
06 Std. Gas	0.081	08:47
07 Room Air	0.000	08:48

Lot No = 01923080A3

Cyl No = 57

Exp Date = 02/05/2025

County = 18

Oper No. = 039947

Remarks:

Annual Inspection Install

Form 106-I8000

Operator Signature JONATHAN LAMPI

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFPD
05/26/2023

Alcohol Analyzer SN 80-006678 8164.14.00 09/16 08:49

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:49
02 Std. Gas	0.081	08:50
03 Room Air	RFI*	08:50
04 Room Air	0.000	08:50

\*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 57

Exp Date = 02/05/2025

County = 18

Remarks:

Oper No. = 039947

Annual Inspection Install

Form 106-I8000

Operator Signature
JONATHAN LAMPI