

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 60 - 60 59 44	CASS Country Jail -1	Machine #1
Reason for Install/Repair	Cas Court 9 Sail	
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	r Location Change	
Dother (Specify) Annual Inspection at	- Crime Laboratory	
Check When Done:	Ŭ	
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Dur	ring Test).
7. Repair and/or Maintenance Performed (if any):	,	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN56		ntoxilyzer® for Use.
3. File Previous Intoxilyzer® Record (SFN504096, Form 120-G		
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date 5/24/23
Crime Laboratory Use Only		,
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crim	
Reviewed/Certified By		Certified Date 20 May 2023

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = CASS 8164.14.00 09/16
05/26/2023 08:07

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? Standard Lot #? Standard Cyl #? 0.080 34921080A1 50 Standard Expiration? 02/05/2024 Oper No? 090801

Flow Cal. Date: 06/20/2013 Slope 677 Intercept -702620

IR Calibration Date: 06/09/2017 3um 9um

Oth Coef(\*100): -16158 -17815
1st Coef(\*100): 266780 136962
2nd Coef(\*100): 1117 1089
H2O adj(mg/l\*10k): 317 364

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature KATHRYN ALLEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944 Location = CASS 8164.14.00 09/16 05/26/2023

08:22

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:22
02 Std. Gas	0.082	08:23
03 Room Air	0.000	08:23
04 Std. Gas	0.081	08:24
05 Room Air	0.000	08:24
06 Std. Gas	0.082	08:25
07 Room Air	0.000	08:25

Lot No = 34921080A1

Cyl No = 50

Exp Date = 02/05/2024

County = 09

Oper No. = 090801

Operator Signature KATHRYN ALLEN

Remarks:
Acceptable

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = CASS 8164.14.00 09/16
05/26/2023 08:26

Test		AC	Time
01 Room Air		0.000	08:28
02 *Subject	Test	RFI*	08:28
03 Room Air		0.000	08:28

\*Invalid Test Inhibited - RFI

Sub Name = TEST INSTAL, TEST INSTAL NA

Sub DOB = 01/01/1900

Sub Sex = Female Weight = 200

Test = DUI Cit = 0200000000000

Dr. Lic. = ND/TEST SUBJECT

Lot No = 34921080A1

Cyl No = 50

Expiration Date = 02/05/2024

County = 09 Oper No. = 090801

operator Signature KATHRYN ALLEN

Remarks:

Form 106-I8000