




**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-004947</b>	Instrument Location <b>Dickinson LEC</b>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Annual Inspection at Crime Laboratory</u>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S). -already done
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): None
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature  <b>Christopher Messer</b>	Date <b>6-1-23</b>
--------------------------------------------------------------------------------------------------------------------------------------------	-----------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By  <b>Anne Wahl</b>	Certified Date <b>02 June 2023</b>
-------------------------------------------------------------------------------------------------------------------------------	---------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004947  
Location = DKPD                8164.14.00 09/16  
06/01/2023                      15:42

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34921080A1  
Standard Cyl #?            40  
Standard Expiration?       02/05/2024  
Oper No?                    010107

Flow Cal. Date:            04/07/2020  
Slope                      618  
Intercept                  -565524

IR Calibration Date:       04/07/2020  
                                    3um                      9um

-----  
0th Coef(\*100):            -22371                      -26151  
1st Coef(\*100):            260822                      133081  
2nd Coef(\*100):            1936                        1144  
H2O adj(mg/l\*10k):        495                         535

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHRISTOPHER MESSER

Remarks:            Print Test

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = DKPD      8164.14.00 09/16  
06/01/2023      15:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:44
02 Std. Gas	0.079	15:45
03 Room Air	0.000	15:45
04 Std. Gas	0.079	15:46
05 Room Air	0.000	15:46
06 Std. Gas	0.078	15:47
07 Room Air	0.000	15:47

Lot No = 34921080A1  
Cyl No = 40  
Exp Date = 02/05/2024  
County = 45      Oper No. = 010107



Operator Signature  
CHRISTOPHER MESSER

Remarks:

ACA Test

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = DKPD      8164.14.00 09/16  
06/01/2023      15:48

Test	AC	Time
01 Room Air	0.000	15:50
02 Reference	RFI*	15:50
03 Room Air	0.000	15:50

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 02/02/2000  
Sub Sex = Male      Weight = 245  
Test = DUI      Cit = 123456  
Dr. Lic. = ND/123456  
Lot No = 34921080A1  
Cyl No = 40  
Expiration Date = 02/05/2024  
County = 45      Oper No. = 010107



Operator Signature  
CHRISTOPHER MESSER

Remarks:

RFI Test

Form 106-I8000