

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004953 Instrument Location Rug	64
Reason for Install/Repair	Ĭ
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change	,
Deter (Specify) Annual Inspection at Crime Laboratory	
0	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	N.
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Durin	ng Test).
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Into	oxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agence	
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Che 104-G), Print Test, ACA Test, and RFI Test.	eckout (SFN59281, Form
Field Inspector Signature	Date /3/7023
Crime Laboratory Use Only	aine alashal concentration from
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determ the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	Eaboratory Division, in the
Reviewed/Certified By	Certified Date 07June2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004953
Location = RUGB	8164.14.00 09/16
06/03/2023	20:41

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
Disable on Memirali:	ies

Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 26021080A1
Standard Cyl #? 8

Standard Cyl #? 8
Standard Expiration? 10/05/2023
Oper No? 130392

Flow Cal. Date: 04/20/2018

Slope 609 Intercept -462839

IR Calibration Date: 04/20/2018
3um 9um

Oth Coef(*100): -22916 -27195
1st Coef(*100): 261768 134580
2nd Coef(*100): 1749 1216
H2O adj(mg/l*10k): 496 540

******* Printer Test End *******

Operator Signature DARREN HEIDBREDER

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = RUGB

Alcohol Analyzer SN 80-004953 8164.14.00 09/16 20:42

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	20:42
02	Std.	Gas	0.080	20:43
03	Room	Air	0.000	20:43
04	Std.	Gas	0.080	20:44
05	Room	Air	0.000	20:44
06	Std.	Gas	0.080	20:45
07	Room	Air	0.000	20:45

Lot No = 26021080A1

Cyl No = 8

06/03/2023

Exp Date = 10/05/2023

County = 35

Oper No. = 130392

Operator Signature DARREN HEIDBREDER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004953
Location = RUGB	8164.14.00 09/16
06/03/2023	20:46

Test	AC	Time
01 Room Air 02 Reference	0.000 RFI*	20:48 20:48
03 Room Air	0.000	20:48

^{*}Invalid Test Inhibited - RFI

Sub Name = TEST RFI, NA NA

Sub DOB = 08/10/1975

Sub Sex = Male

Weight = 120

Test = DUI

Cit = 1

Dr. Lic. = ND/12345678

Lot No = 26021080A1

Cyl No = 8

Expiration Date = 10/05/2023

County = 35

Oper No. = 130392

Operator Signature DARREN HEIDBREDER

Remarks:

Form 106-I8000