

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 6509	Instrument Location Grand Forks AFB ND					
Reason for Install/Repair						
☐ Install After Receiving From Crime Laboratory ☐ Install Aft	er Location Change					
Dother (Specify) Annual Inspection a	+ Crime laboratory					
	0					
Check When Done:						
1. Surge Protector Installed/Property Grounded.						
Z Telephone Line Connected to Intoxilyzer® 8000.						
3. Breath Tube Heated.						
Enter Preliminary Data (i.e. Date, Time, DST (Y), and Local	tion; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).						
6. Run Tests:						
A. Print Test (Level1, Function P).						
B. ACA Test (Level 1, Function C).						
C. Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):						
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form					
Field Inspector Siggature	Date .					
O miles	06 Jun 23					
Scient to the Sci						
Crime Laboratory Use Only						
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct of	be used for the analysis of breath to determine alcohol concentration from file at the Office of Attorney General, Crime Laboratory Division, in the opy of the documents received.					
Reviewed/Certified By	Certified Date 08. June 2022					
	0 0001 42025					

From Maran 14,2023 Chain of custody, no install papernonk was received. Upon connecting to the instrument, it was found that no CMS test records were ran.

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006509
Location = GFAF 8164.14.00 09/16
06/06/2023 13:38

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Flow Cal. Date: 06/11/2019 Slope 445 Intercept 431160

IR Calibration Date: 06/11/2019
3um 9um

Oth	Coef(*100):	-27266	-14145
1st	Coef(*100):	276865	139368
2nd	Coef(*100):	2326	1158
H20	adj $(mg/l*10k)$:	584	288

****** Printer Test End *******

Operator Signature MATTHEW OXENFIELD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFAF
06/06/2023

Alcohol Analyzer SN 80-006509 8164.14.00 09/16 13:39

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	13:39
02	Std.	Gas	0.079	13:40
03	Room	Air	0.000	13:40
04	Std.	Gas	0.080	13:41
05	Room	Air	0.000	13:41
06	Std.	Gas	0.079	13:41
07	Room	Air	0.000	13:42

Lot No = 01923080A3

Cyl No = 58

Exp Date = 02/05/2025

County = 18

Oper No. = 133391

Operator Signature MATTHEW OXENFIELD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFAF
06/06/2023

Alcohol Analyzer SN 80-006509 8164.14.00 09/16 13:43

DRY CAL CHECK

Te	est		AC	Time
01	Room	Air	RFI*	13:44
02	Room	Air	0.000	13:45

*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 58

Exp Date = 02/05/2025

County = 18

Oper No. = 133391

Operator Signature MATTHEW OXENFIELD

Remarks:

Form 106-I8000