

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	damstairs			
Serial Number 80 - 00 650 2	Instrument Location STUTSMAN CO Corr = STUT(JUV)			
Reason for Install/Repair				
Install After Receiving From Crime Laboratory				
Dether (specify) Annual Inspection at Chine Laboratory				
Check When Done:	0			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
🔄 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature	Date 6-7-23			

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	~		
Reviewed/Certified By	Am	Mall	Certified Date 08June2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006502 Location = STUT 8164.14.00 09/16 06/07/2023 22:30 ************ Printer Test **************** abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 26021080A1 Standard Cyl #? 5 Standard Expiration? 10/05/2023 Oper No? 132539 Flow Cal. Date: 08/19/2015 Slope 674 Intercept -654491 IR Calibration Date: 08/19/2015 3um 9um _____ 0th Coef(*100): -24446 -181031st Coef(*100): 267787 2nd Coef(*100): 3387 1st Coef(*100): 135417 1610 H2O adj(mg/l*10k): 740 513 ********** Printer Test End ***********

Operator Signature JOSHUA FASTNACHT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-006502Location = STUT
06/07/20238164.14.00 09/16
22:32

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.081 0.000 0.080 0.000	22:32 22:33 22:33 22:34 22:34 22:35 22:35

Lot No = 26021080A1 Cyl No = 5 Exp Date = 10/05/2023 County = 47

Oper No. = 132539

D Operator Signature

JOSHUA FASTNACHT

Remarks:

Form 106-18000

•

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-006502Location = STUT
06/07/20238164.14.00 09/16
22:36

DRY CA	L CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* 0.000	22:36 22:37
*Invalid Test Inhibited - RFI		
Lot No = 26021080A1 Cyl No = 5 Exp Date = 10/05/2023	3	

Oper No. = 132539

 \sim

Operator Signature JOSHUA FASTNACHT

Remarks:

County = 47

Form 106-18000

r