

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	5 5
Reason for Install/Repair	Grafton - Walsh	Co Sherilf's
	er Location Change	
Dother (Specify) Annual Inspection at	- Crime Laboratory	
	0	
Check When Done:		
Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode o	r Level 1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Place it by the I	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair 0	Checkout (SFN59281, Form
Field Inspector Signature	<u>-</u>	Date
471	\$	6-15-93
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to	be used for the analysis of breath to dete	ermine alcohol concentration from
the date the Field Inspector performed the installation. This record on	file at the Office of Attorney General, Cri	
County of Burleigh, North Dakota, is certified to be a true and correct of	opy of the documents received.	
Reviewed/Certified By		Certified Date
		201410 12022

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004948
Location = GRAF	8164.14.00 09/16
06/12/2023	15:38

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFĞHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current	Instrument	Setup
	try Mode:	-

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	01923080A3
Standard Cyl #?	21
Standard Expiration?	02/05/2025
Oper No?	132512

Flow Cal. Date:	08/28/2018
Slope	642
Intercept	-606282

IR Calibration Date	e: 08/28/2	018
	3um	9um
0th Coef(*100):	-13609	-14548
1st Coef(*100):	248186	128726
2nd Coef(*100):	3538	1508
H20 adj(mg/l*10k):	602	515

******* Printer Test End *******

Operator Signature JORDAN HOFFER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GRAF
06/15/2023

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 09:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:37
02 Std. Gas	0.082	09:37
03 Room Air	0.000	09:38
04 Std. Gas	0.082	09:38
05 Room Air	0.000	09:39
06 Std. Gas	0.082	09:39
07 Room Air	0.000	09:40

Lot No = 01923080A3

Cyl No = 21

Exp Date = 02/05/2025

County = 50

Oper No. = 132512

Operator Signature JORDAN HOFFER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004948
Location = GRAF	8164.14.00 09/16
06/15/2023	09:40

Test	AC	Time
01 Room Air	0.000	09:42
02 *Subject Test	RFI*	09:42
03 Room Air	0.000	09:43

*Invalid Test Inhibited - RFI

Sub Name = J, J K Sub DOB = 11/11/1990Sub Sex = Male Weight = 456Test = DUI Cit = NA Dr. Lic. = ND/NALot No = 01923080A3Cyl No = 21Expiration Date = 02/05/2025County = 50

Oper No. = 132512

Operator Signature JORDAN HOFFER

Remarks:

Form 106-I8000