

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Coriel Number				
Serial Number 80-00667/	Instrument Locations Lake Region Cor	r = DEVL		
Reason for Install/Repair				
Install After Receiving From Crime Laboratory Install After	er Location Change			
Decity) Annual Inspection at	- Crime Laboratory_			
Check When Done:	Ö			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	,			
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).		
7. Repair and/or Maintenance Performed (if any): Ahnag(1 tracetion			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the I	ntoxilvzer® for Use		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date		
A Marine May	539	6/17/2023		
Crime Laboratory Use Only	,			
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin	rmine alcohol concentration from ne Laboratory Division, in the		
Reviewed/Certified By Mul		Certified Date 27June2023		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006671
Location = DEVL 8164.14.00 09/16
06/17/2023 20:08

********* Printer Test *********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 19 Standard Expiration? 02/05/2025 Oper No? 132008

Flow Cal. Date: 08/22/2016 Slope 703 Intercept -691868

IR Calibration Date: 08/22/2016
3um 9um

Oth Coef(*100): -17135 -24345
1st Coef(*100): 254798 132856
2nd Coef(*100): 3386 1354
H2O adj(mg/l*10k): 610 539

******* Printer Test End ********

Operator Signature

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL
06/17/2023

SN 80-006671 8164.14.00 09/16 20:01

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	20:02
02	Std.	Gas	0.080	20:02
03	Room	Air	0.000	20:03
04	Std.	Gas	0.080	20:03
05	Room	Air	0.000	20:03
06	Std.	Gas	0.080	20:04
07	Room	Air	0.000	20:04

Lot No = 01923080A3

Cyl No = 19

Exp Date = 02/05/2025

County = 36

Oper No. = 132008

Operator Signature
WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL
06/17/2023

Alcohol Analyzer SN 80-006671 8164.14.00 09/16 20:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:05
02 Std. Gas	RFI*	20:06
03 Room Air	0.000	20:06

*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 19

Exp Date = 02/05/2025

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000