

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

T	T. T			
Serial Number 80 - 00 30 6 7	Instrument Location Grafton	1 Pd		
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change			
Dother (Specify) Annual Inspection at	- Crime Laboratory			
	0			
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.		
9 File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature  Date  Le[19/2023				
Crime Laboratory Use Only		, , =		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin			
Reviewed/Certified By		Certified Date 20June2023		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = GFTN 8164.14.00 09/16
06/19/2023 09:47

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	34921080A1
Standard Cyl #?	36
Standard Expiration?	02/05/2024
Oper No?	020305

Flow Cal. Date:	04/17/2020
Slope	609
Intercept	-236750

IR	Calibration	Date:	04/17/2020
		3um	9um

0th	Coef(*100):	-19307	-18350
	Coef(*100):	282277	142540
2nd	Coef(*100):	1257	744
H20	adj(mg/l*10k):	302	214

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
ANTHONY DUMAS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
06/19/2023

Alcohol Analyzer SN 80-003067 8164.14.00 09/16 09:57

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:58
02 Std. Gas	0.077	09:58
03 Room Air	0.000	09:59
04 Std. Gas	0.077	09:59
05 Room Air	0.000	10:00
06 Std. Gas	0.078	10:00
07 Room Air		10:01

Lot No = 34921080A1

Cyl No = 36

Exp Date = 02/05/2024

County = 50

Oper No. = 020305

Operator Signature ANTHONY DUMAS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
06/19/2023

SN 80-003067 8164.14.00 09/16 10:05

Test			AC	Time	
01	Room	Air		RFI*	10:06
02	Room	Air	2.5	RFI*	10:06

\*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T Sub DOB = 01/01/1900

Sub Sex = Male

Weight = 200 Cit = 1234567

Test = OTH
Dr. Lic. = ND/TES001234

Lot No = 34921080A1

Cyl No = 36

Expiration Date = 02/05/2024

County = 50

Oper No. = 020305

Operator Signature ANTHONY DUMAS

Remarks:

RFI Test after annual Insp.

Form 106-I8000