



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION**

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-00 5357 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No

B. Calibration Adjustment (Level 3,M,C,O)

1. Autocalibration Printout Attached
 - Max Power Res Value ≥ 10
 - Auto Range Res Value ≥ 4
2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H ₂ O	NA-Milli-Q H ₂ O	MP3066
2	0.040	202111A	09Nov23	MP6040
3	0.080	202110C	26Oct23	MP5320
4	0.100	202304A	04Apr25	MP5290
5	0.300	202201F	18Jan24	MP3059

3. 0.080 AC Calibration Gas for H₂O Adjustment

Lot No. 26021080A Cyl No. 20 Exp. Date: 10/5/23

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 961 mbar
 Adjusted to using barometer 958 mbar
 Auto Calibration Report printout 961 mbar
 Barometer Model 03316-72
 Barometer Serial Number 881001
 Barometer Calibration Expiration Date 9/1/23

5. Screen displayed "Calibration Success"
6. Calibration Adjustment Printout Attached
 - Solution 1 Avg % Abs ≤ 0.2500
 - Solution 2-5 REL STD DEV ≤ 3.000

Residual (g/210 L) values for solutions 1 - 5 ≤ 0.0020 for 3 μm and 9 μm channels

Dry Gas H₂O adjustment sum for 3 μm and 9 μm channels within ± 10

3 μm 3494 (Ave.) + 315 (H₂O Adj.) = 3809

9 μm 3391.3 (Ave.) + 418 (H₂O Adj.) = 3809.3

C. Is an Annual Inspection due for this instrument? Yes or No

If Yes, complete Intoxilyzer 8000 Annual Inspection (Qualtrax ID: 11698)

If No, complete Intoxilyzer 8000 Calibration (Qualtrax ID: 11871).

Remarks/Notes: N/A


Analyst Signature

21 June 2023
Date


Reviewer Signature

22 June 2023
Date

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005357
06/21/2023 11:24:26

Auto Calibration
Max Power Res Value = 57
Auto Range Res Value = 38

Auto Calibration Printout
An Hall

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005357
 06/21/2023 11:24:26

Auto Calibration

pg 1 of 2

```

<<<<<      3um      >>>>>      <<<<<      9um      >>>>>
-----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     0.0610      (0.0100)          0.1160      (-0.0060)
Sample #2     0.0250      (0.0630)          0.0760      (0.0310)
Sample #3     0.0100      (0.0940)          0.1100      (0.0250)
Sample #4     0.0540      (0.1000)          0.1170      (0.0240)
Avg % Abs     0.0297      (0.0857)          0.1010      (0.0267)
STD DEV       0.0224      (0.0199)          0.0219      (0.0038)
REL STD DEV   75.398      (23.180)          21.715      (14.197)
  
```

```

-----
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     0.7450      (-0.0110)         1.4030      (-0.0180)
Sample #2     0.7440      (0.0300)          1.4830      (-0.0330)
Sample #3     0.7480      (0.0100)          1.4580      (-0.0310)
Sample #4     0.7380      (0.0360)          1.4640      (-0.0230)
Avg % Abs     0.7433      (0.0253)          1.4683      (-0.0290)
STD DEV       0.0050      (0.0136)          0.0131      (0.0053)
REL STD DEV   0.677       (53.738)          0.889       (18.247)
  
```

```


-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     1.4350      (0.0190)          2.7580      (0.0200)
Sample #2     1.4660      (0.0120)          2.7660      (0.0530)
Sample #3     1.4940      (0.0070)          2.8070      (0.0420)
Sample #4     1.4810      (0.0120)          2.8030      (0.0380)
Avg % Abs     1.4803      (0.0103)          2.7920      (0.0443)
STD DEV       0.0140      (0.0029)          0.0226      (0.0078)
REL STD DEV   0.947       (27.936)          0.810       (17.521)
  
```

```

-----
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     1.7800      (0.0060)          3.4400      (0.0090)
Sample #2     1.8080      (0.0180)          3.4840      (0.0500)
Sample #3     1.8470      (0.0040)          3.5030      (0.0310)
Sample #4     1.8480      (0.0010)          3.4590      (0.0450)
Avg % Abs     1.8343      (0.0077)          3.4820      (0.0420)
STD DEV       0.0228      (0.0091)          0.0221      (0.0098)
REL STD DEV   1.244       (118.354)         0.634       (23.450)
  
```

```

-----
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1     5.1300      (0.0000)          9.4220      (-0.0010)
Sample #2     5.2370      (0.0240)          9.6050      (0.0360)
Sample #3     5.2440      (0.0200)          9.5970      (0.0490)
Sample #4     5.2170      (0.0150)          9.5640      (0.0460)
Avg % Abs     5.2327      (0.0197)          9.5887      (0.0437)
STD DEV       0.0140      (0.0045)          0.0217      (0.0068)
REL STD DEV   0.268       (22.928)          0.227       (15.588)
  
```



TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005357
 06/21/2023 11:24:26

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

 Zero Order Coef -60.39
 First Order Coef 2574.61
 Second Order Coef 31.87

<<<<< 9um >>>>>

 Zero Order Coef -130.97
 First Order Coef 1360.71
 Second Order Coef 14.88

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0003
0.040	0.039	0.0007
0.080	0.080	-0.0002
0.100	0.100	-0.0002
0.300	0.300	0.0000

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0001
0.040	0.040	0.0001
0.080	0.079	0.0005
0.100	0.101	-0.0005
0.300	0.300	0.0000

<<<<< 3um >>>>>

<<<<< 9um >>>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

Sample	3um	9um
Sample #1	3565.00	3428.00
Sample #2	3520.00	3386.00
Sample #3	3420.00	3377.00
Sample #4	3542.00	3411.00
Avg	3494.0000	3391.3333
STD DEV	65.0231	17.6163
REL STD DEV	1.861	0.519
H2O adjust (mg/l*10k)	315	418

Atmospheric Pressure = 961

*****CALIBRATION SUCCESSFUL*****

