



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION**

**INTOXILYZER® 8000 CALIBRATION ADJUSTMENT**

Intoxilyzer® 8000 Serial Number: 80-00 4206 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No

B. Calibration Adjustment (Level 3,M,C,O)

1.  Autocalibration Printout Attached
  - Max Power Res Value  $\geq 10$
  - Auto Range Res Value  $\geq 4$
2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H <sub>2</sub> O	NA-Milli-Q H <sub>2</sub> O	MP3066
2	0.040	202111A	09NOV23	MP6040
3	0.080	202110C	26 Oct 23	MP5320
4	0.100	202304A	04 Apr 25	MP5290
5	0.300	202201F	18 Jan 24	MP3059

3. 0.080 AC Calibration Gas for H<sub>2</sub>O Adjustment

Lot No. 26021080A Cyl No. 20 Exp. Date: 10/5/23

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 949 mbar  
 Adjusted to using barometer 963 mbar  
 Auto Calibration Report printout 963 mbar  
 Barometer Model 03316-72  
 Barometer Serial Number 881001  
 Barometer Calibration Expiration Date 9/1/23

5.  Screen displayed "Calibration Success"
6.  Calibration Adjustment Printout Attached
  - Solution 1 Avg % Abs  $\leq 0.2500$
  - Solution 2-5 REL STD DEV  $\leq 3.000$

Residual (g/210 L) values for solutions 1 - 5  $\leq 0.0020$  for 3  $\mu$ m and 9  $\mu$ m channels

Dry Gas H<sub>2</sub>O adjustment sum for 3  $\mu$ m and 9  $\mu$ m channels within  $\pm 10$

3  $\mu$ m 3338.6 (Ave.) + 471 (H<sub>2</sub>O Adj.) = 3809.6

9  $\mu$ m 3189.6 (Ave.) + 620 (H<sub>2</sub>O Adj.) = 3809.6

C. Is an Annual Inspection due for this instrument? Yes or No

If Yes, complete Intoxilyzer 8000 Annual Inspection (Qualtrax ID: 11698)


If No, complete Intoxilyzer 8000 Calibration (Qualtrax ID: 11871).

Remarks/Notes: NIA

---

---

---

  
Analyst Signature

21 June 2023  
Date

  
Reviewer Signature

22 June 2023  
Date

TOXL  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-004206  
06/21/2023 13:18:58

Auto Calibration  
Max Power Res Value = 24  
Auto Range Res Value = 8

Auto Calibration printout

*Oh*

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004206  
 06/21/2023 13:18:58

Auto Calibration

pg 1 of 2


	<<<< 3um >>>>		<<<< 9um >>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1				
Sample #1	0.1090	(-0.0290)	0.1540	(-0.0020)
Sample #2	0.0450	(0.0590)	0.1550	(0.0330)
Sample #3	0.0880	(0.0880)	0.1660	(0.0300)
Sample #4	0.0510	(0.1330)	0.1790	(0.0560)
Avg % Abs	0.0613	(0.0933)	0.1667	(0.0397)
STD DEV	0.0233	(0.0373)	0.0120	(0.0142)
REL STD DEV	37.970	(39.951)	7.208	(35.860)

	<<<< 3um >>>>		<<<< 9um >>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1				
Sample #1	0.7210	(0.0100)	1.5630	(0.0050)
Sample #2	0.7670	(0.0070)	1.5660	(0.0270)
Sample #3	0.7620	(0.0340)	1.5900	(0.0110)
Sample #4	0.7680	(0.0440)	1.5940	(0.0230)
Avg % Abs	0.7657	(0.0283)	1.5833	(0.0203)
STD DEV	0.0032	(0.0191)	0.0151	(0.0083)
REL STD DEV	0.420	(67.552)	0.956	(40.951)

	<<<< 3um >>>>		<<<< 9um >>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1				
Sample #1	1.4370	(0.0110)	2.9070	(0.0080)
Sample #2	1.4810	(0.0280)	2.9930	(0.0210)
Sample #3	1.4750	(0.0370)	2.9710	(0.0370)
Sample #4	1.4910	(0.0290)	2.9810	(0.0350)
Avg % Abs	1.4823	(0.0313)	2.9817	(0.0310)
STD DEV	0.0081	(0.0049)	0.0110	(0.0087)
REL STD DEV	0.545	(15.743)	0.369	(28.122)

	<<<< 3um >>>>		<<<< 9um >>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1				
Sample #1	1.8030	(0.0160)	3.6040	(0.0100)
Sample #2	1.8350	(0.0160)	3.7220	(0.0040)
Sample #3	1.8630	(0.0300)	3.7040	(0.0330)
Sample #4	1.8340	(0.0250)	3.6670	(0.0370)
Avg % Abs	1.8440	(0.0237)	3.6977	(0.0247)
STD DEV	0.0165	(0.0071)	0.0280	(0.0180)
REL STD DEV	0.893	(29.977)	0.758	(73.010)

	<<<< 3um >>>>		<<<< 9um >>>>	
	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1				
Sample #1	5.0840	(0.0190)	9.8720	(0.0250)
Sample #2	5.2560	(0.0080)	10.1770	(0.0300)
Sample #3	5.3080	(0.0130)	10.2100	(0.0510)
Sample #4	5.2400	(0.0320)	10.1660	(0.0650)
Avg % Abs	5.2680	(0.0177)	10.1843	(0.0487)
STD DEV	0.0356	(0.0127)	0.0229	(0.0176)
REL STD DEV	0.675	(71.673)	0.225	(36.198)



AEN

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004206  
 06/21/2023 13:18:58

Auto Calibration

pg 2 of 2

<<<< 3um >>>>			<<<< 9um >>>>		
-----			-----		
Zero Order Coef	-148.68			-209.07	
First Order Coef	2636.67			1308.72	
Second Order Coef	19.57			11.23	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	0.000	-0.0003	0.000	0.000	-0.0002
0.040	0.040	0.0005	0.040	0.040	0.0003
0.080	0.080	0.0001	0.080	0.080	0.0003
0.100	0.100	-0.0004	0.100	0.100	-0.0005
0.300	0.300	0.0000	0.300	0.300	0.0000
-----			-----		

<<<< 3um >>>>			<<<< 9um >>>>		
-----			-----		
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1					
Sample					
Sample #1		3378.00		3148.00	
Sample #2		3383.00		3211.00	
Sample #3		3243.00		3153.00	
Sample #4		3390.00		3205.00	
Avg		3338.6667		3189.6667	
STD DEV		82.9237		31.8957	
REL STD DEV		2.484		1.000	
H2O adjust (mg/l*10k)		471		620	

Atmospheric Pressure = 963

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*