



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION**

**INTOXILYZER® 8000 CALIBRATION ADJUSTMENT**

Intoxilyzer® 8000 Serial Number: 80-00 5956 Calibration Adjustment Location: TOXL

**A. Pre-Adjustment**

Replaced Simulator Return O-Ring Yes or No

**B. Calibration Adjustment (Level 3,M,C,O)**

1.  Autocalibration Printout Attached
  - Max Power Res Value  $\geq 10$
  - Auto Range Res Value  $\geq 4$
2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H <sub>2</sub> O	NA-Milli-Q H <sub>2</sub> O	MP3046
2	0.040	2022111A	09 Nov 23	MP6040
3	0.080	202110C	26 Oct 23	MP5320
4	0.100	202304A	04 Apr 25	MP5290
5	0.300	202201F	18 Jan 24	MP3059

3. 0.080 AC Calibration Gas for H<sub>2</sub>O Adjustment

Lot No. 26021080A Cyl No. 20 Exp. Date: 10/5/23

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 955 mbar  
 Adjusted to using barometer 961 mbar  
 Auto Calibration Report printout 961 mbar  
 Barometer Model 02316-72  
 Barometer Serial Number 881001  
 Barometer Calibration Expiration Date 9/1/23

5.  Screen displayed "Calibration Success"
6.  Calibration Adjustment Printout Attached
  - Solution 1 Avg % Abs  $\leq 0.2500$
  - Solution 2-5 REL STD DEV  $\leq 3.000$

TOXL  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-005956  
06/21/2023 13:57:48

Auto Calibration  
Max Power Res Value = 24  
Auto Range Res Value = 8

Auto Calibration Printout

*A. Smith*

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005956  
 06/21/2023 13:57:48

Auto Calibration

pg 1 of 2

```

<<<<<      3um      >>>>>
-----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1 0.0310      (0.0160)      0.1010      (0.0080)
Sample #2 0.0610      (0.0550)      0.1220      (0.0050)
Sample #3 0.0430      (0.0920)      0.1350      (0.0260)
Sample #4 0.0640      (0.1030)      0.1390      (0.0330)
Avg % Abs 0.0560      (0.0833)      0.1320      (0.0213)
STD DEV   0.0114      (0.0251)      0.0089      (0.0146)
REL STD DEV 20.282      (30.175)      6.733       (68.305)
  
```

```

-----
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1 0.6910      (-0.0020)      1.4270      (0.0040)
Sample #2 0.7040      (0.0230)      1.4470      (0.0370)
Sample #3 0.6970      (0.0380)      1.4820      (0.0200)
Sample #4 0.7070      (0.0300)      1.5150      (0.0050)
Avg % Abs 0.7027      (0.0303)      1.4813      (0.0207)
STD DEV   0.0051      (0.0075)      0.0340      (0.0160)
REL STD DEV 0.730       (24.744)      2.296       (77.470)
  
```

```


-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1 1.3460      (0.0200)      2.7460      (0.0160)
Sample #2 1.3600      (0.0480)      2.8060      (0.0290)
Sample #3 1.4080      (0.0320)      2.8540      (0.0230)
Sample #4 1.3850      (0.0530)      2.8490      (0.0330)
Avg % Abs 1.3843      (0.0443)      2.8363      (0.0283)
STD DEV   0.0240      (0.0110)      0.0264      (0.0050)
REL STD DEV 1.734       (24.744)      0.930       (17.764)
  
```

```

-----
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1 1.7060      (-0.0190)      3.3990      (-0.0280)
Sample #2 1.7500      (-0.0150)      3.4990      (-0.0050)
Sample #3 1.7100      (0.0170)      3.5300      (-0.0170)
Sample #4 1.7170      (0.0040)      3.5120      (-0.0230)
Avg % Abs 1.7257      (0.0020)      3.5137      (-0.0150)
STD DEV   0.0214      (0.0161)      0.0156      (0.0092)
REL STD DEV 1.238       (804.674)      0.443       (61.101)
  
```

```

-----
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1
Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1 4.7550      (-0.0140)      9.3630      (-0.0110)
Sample #2 4.9020      (0.0090)      9.6360      (0.0050)
Sample #3 4.9020      (0.0310)      9.6920      (0.0380)
Sample #4 4.9020      (0.0390)      9.6900      (0.0300)
Avg % Abs 4.9020      (0.0263)      9.6727      (0.0243)
STD DEV   0.0000      (0.0155)      0.0318      (0.0172)
REL STD DEV 0.000       (58.993)      0.328       (70.744)
  
```



AEN

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-005956  
 06/21/2023 13:57:48

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -129.62
First Order Coef  2803.78
Second Order Coef  27.83
  
```

```

<<<<<      9um      >>>>>
-----
Zero Order Coef   -165.16
First Order Coef  1359.18
Second Order Coef  13.92
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.001    -0.0006
0.040    0.039    0.0011
0.080    0.080    0.0001
0.100    0.101    -0.0006
0.300    0.300    0.0001
-----
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.000    -0.0003
0.040    0.039    0.0005
0.080    0.080    0.0002
0.100    0.100    -0.0004
0.300    0.300    0.0000
-----
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
  
```

```

Sample #1          3329.00          3307.00
Sample #2          3360.00          3320.00
Sample #3          3453.00          3344.00
Sample #4          3310.00          3335.00
Avg                3374.3333          3333.0000
STD DEV            72.5695           12.1244
REL STD DEV        2.151             0.364
H2O adjust (mg/l*10k) 435             476
  
```

Atmospheric Pressure = 961

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*