

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 005955	Instrument Location UNDPD			
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change				
Dother (Specify) Annual Inspection at	- Crime Laboratory			
Check When Done:	0			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
RIM 6. Run Tests:				
TM MA. Print Test (Level1, Function P).				
TM B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
ルM 図 7. Repair and/or Maintenance Performed (if any): None	8.			
72M 38. Complete the Top Portion of the Intoxilyzer® Record (SFN50				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature		Date 6-30-23		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By		Certified Date 03Jul2023		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005955
Location = UNDP 8164.14.00 09/16
06/30/2023 12:38

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Prellm Asic.
Display Third Digit?
Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull?
of Print Copies? Yes 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34921080A1 30 Standard Expiration? 02/05/2024 119409 Oper No?

Flow Cal. Date: 04/19/2022 Slope 2446 Intercept -3112812

IR Calibration Date: 06/21/2023 3um 9um

Oth Coef(*100): -12629 -18041
1st Coef(*100): 273918 134302
2nd Coef(*100): 3542 1308
H2O adj(mg/l*10k): 397 389

****** Printer Test End *******

perator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005955 Location = UNDP 06/30/2023

8164.14.00 09/16 12:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:47
02 Std. Gas	0.079	12:48
03 Room Air	0.000	12:48
04 Std. Gas	0.078	12:48
05 Room Air	0.000	12:49
06 Std. Gas	0.079	12:49
07 Room Air	0.000	12:50

Lot No = 34921080A1

Cyl No = 30

Exp Date = 02/05/2024

County = 18

Oper No. = 119409

Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005955 Location = UNDP 8164.14.00 09/16 06/30/2023

Test	AC	Time
01 Room Air	0.000	12:53
02 *Subject Test	RFI*	12:53
03 Room Air	0.000	12:54

12:52

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 01/02/1967

Weight = 185 Sub Sex = MaleTest = DUI Cit = 00000000

Dr. Lic. = ND/0000000

Lot No = 34921080A1

Cyl No = 30

Expiration Date = 02/05/2024

Oper No. = 119409County = 18

> Operator Signature TRACY MEIDINGER

Remarks:

Form 106-I8000