

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number SO - 006501 Instrument Location Mino+ AFI3				
Reason for Install/Repair				
Install After Receiving From Crime Laboratory   Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
★ 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).			
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the In	ntoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ency.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form			
Field Inspector Signature	Date			
	13 July 2023			
Crime Laboratory Use Only	Ü			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the				
County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By	Cartified Date			
	Certified Date			

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006501
Location = MNAF 8164.14.00 09/16
07/13/2023 15:00

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34921080A1 43 Standard Expiration? 02/05/2024 Oper No? 133426

Flow Cal. Date: 06/13/2019 Slope 651 Intercept -777821

IR Calibration Date: 06/13/2019
3um 9um

Oth Coef(\*100): -20531 -14906 1st Coef(\*100): 270252 139749 2nd Coef(\*100): 1917 1284 H2O adj(mg/l\*10k): 393 321

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CODY BAKER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = MNAF 8164.14.00 09/16 07/13/2023

Alcohol Analyzer SN 80-006501 15:01

Test			AC	Time
01	Room	Air	RFI*	15:05
02	Room	Air	RFI*	15:05

\*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST NA

Sub DOB = 01/01/1999

Sub Sex = Unknown Weight = NA Test = OTH Cit = NA

Dr. Lic. = ND/NA

Lot No = 34921080A1

Cyl No = 43

Expiration Date = 02/05/2024

County = 64Oper No. = 133426

Operator Signature CODY BAKER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MNAF
06/23/2023

SN 80-006501 8164.14.00 09/16 09:47

## DRY CAL CHECK

Test	AC		Time
01 Room Air	0.000		09:48
02 Std. Gas	0.078		09:48
03 Room Air	0.000		09:49
04 Std. Gas	0.078	1	09:49
05 Room Air	0.000		09:50
06 Std. Gas	0.077		09:50
07 Room Air	0.000		09:50

Lot No = 34921080A1

Cyl No = 43

Exp Date = 02/05/2024

County = 64

Oper No. = 133290

Operator Signature BRYCE SLOAN

Remarks:

Form 106-I8000