

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-001096 Instrument Location Beau And Police De	Dagiment		
Reason for Install/Repair	FORT ITSUIT		
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change			
Decity Annual Inspection at Crime Laboratory			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).		
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Ir	itoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.		
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form		
Field Inspector Signature WWW.742	Date 5-12-2023		
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.			
Reviewed/Certified By	Certified Date 15Aug2023		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-007096
Location = BEPD	8164.16.00 09/18
05/12/2023	15:20

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
CI 1 1 TT- 1 0	0 000

Standard Value? 0.080
Standard Lot #? 26021080A1
Standard Cyl #? 36
Standard Expiration? 10/05/2023
Oper No? 130634

10/02/2018 Flow Cal. Date: 695 Slope -692605 Intercept

IR Calibration Date: 10/02/2018
3um 9um

_____ Oth Coef(*100): -33544 -18475 1st Coef(*100): 261464 138054 2nd Coef(*100): 3020 1612 H20 adj(mq/l*10k): 573

******* Printer Test End *******

Operator Signature WILLIAM ELMORE

WKlst742

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BEPD
05/12/2023

Alcohol Analyzer SN 80-007096 8164.16.00 09/18 15:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:22
02 Std. Gas	0.079	15:22
03 Room Air	0.000	15:23
04 Std. Gas	0.080	15:23
05 Room Air	0.000	15:24
06 Std. Gas	0.079	15:24
07 Room Air	0.000	15:24

Lot No = 26021080A1

Cyl No = 36

Exp Date = 10/05/2023

County = 29

Oper No. = 130634

Operator Signature
WILLIAM ELMORE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BEPD
05/12/2023

Alcohol Analyzer SN 80-007096 8164.16.00 09/18 15:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	15:26
02 Room Air	RFI*	15:26

*Invalid Test Inhibited - RFI

Lot No = 26021080A1

Cyl No = 36

Exp Date = 10/05/2023

County = 29

Oper No. = 130634

Operator Signature WILLIAM ELMORE

Remarks:

Form 106-I8000