

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005943 Instrument Location NCVILL hospital						
Reason for Install/Repair						
Install After Receiving From Crime Laboratory						
Other (Specify)						
Check When Done:						
Surge Protector Installed/Property Grounded.						
2. Telephone Line Connected to Intoxilyzer® 8000.						
🖄 3. Breath Tube Heated.						
🖄 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).						
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).						
6. Run Tests:						
A. Print Test (Level1, Function P).						
B. ACA Test (Level 1, Function C).						
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).						
7. Repair and/or Maintenance Performed (if any):						
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.						
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.						
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.						
Field Inspector Signature Date						

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005943 Location = NCSO 8164.14.00 09/16 08/04/2023 13:45 *********** Printer Test ************ abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Dry 0.080 Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? 34921080A1 26 Standard Expiration? 02/05/2024 Oper No? 100754 06/20/2013 Flow Cal. Date: Slope 686 -735246 Intercept IR Calibration Date: 08/26/2016 3um 9um _____ Oth Coef(*100): -13589 -19547
 1st Coef(*100):
 284566

 2nd Coef(*100):
 2717
139614 1336 H2O adj(mg/l*10k): 381 404

Operator Signature

DANIELLE BJORLIE

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-005943Location = NCSO8164.14.00 09/16
13:45

			DRY CAL CHECK	
Τe	est		AC	Time
02 03 04 05	Room Std. Room Std. Room Std. Room	Gas Air Gas Air Gas	0.000 0.079 0.000 0.080 0.000 0.000 0.079 0.000	13:46 13:47 13:47 13:47 13:48 13:48 13:48 13:49

Lot No = 34921080A1 Cyl No = 26 Exp Date = 02/05/2024 County = 32 Oper No. = 100754

Operator Signature DANIELLE BJORLIE

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-005943Location = NCSO8164.14.00 09/16
13:49

DRY CAL		
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air *Invalid Test Inhibited - RFI	0.000 RFI* 0.000	13:50 13:50 13:51
Lot No = $34921080A1$ Cyl No = 26 Exp Date = $02/05/2024$ County = 32	Oper No	= 100754

Operator Signature DANIELLE BJORLIE

Remarks:

Form 106-18000