



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59281 (06/2018)

Serial Number <b>80-004199</b>	Instrument Location <b>Cass County Jail-Machine #2</b>
Reason for Install/Repair	
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

### Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <b>[Signature]</b>	Date <b>8-26-23</b>
---	------------------------

### Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <b>[Signature]</b>	Certified Date <b>29Aug2023</b>
---	------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004199  
Location = CASS                    8164.14.00 09/16  
08/26/2023                            00:03

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            01923080A3  
Standard Cyl #?            23  
Standard Expiration?        02/05/2025  
Oper No?                    090801

Flow Cal. Date:            06/17/2019  
Slope                        593  
Intercept                   -481610

IR Calibration Date:        06/17/2019  
                                  3um                    9um

-----  
0th Coef(\*100):            -17512                -19237  
1st Coef(\*100):            277207                136841  
2nd Coef(\*100):            1620                   1027  
H2O adj(mg/l\*10k):        211                    375

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
KATHRYN ALLEN

Remarks:

*Reinstal  
at Tail !*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = CASS      8164.14.00 09/16  
08/26/2023      00:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	00:04
02 Std. Gas	0.082	00:04
03 Room Air	0.000	00:05
04 Std. Gas	0.081	00:05
05 Room Air	0.000	00:06
06 Std. Gas	0.082	00:06
07 Room Air	0.000	00:07

Lot No = 01923080A3  
Cyl No = 23  
Exp Date = 02/05/2025  
County = 09      Oper No. = 090801

  
Operator Signature  
KATHRYN ALLEN

Remarks:

Form 106-I8000

*Instal  
back at  
Jan.*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = CASS      8164.14.00 09/16  
08/26/2023      00:08

Test	AC	Time
01 Room Air	0.000	00:09
02 *Subject Test	RFI*	00:09
03 Room Air	0.000	00:09

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 01/01/1900  
Sub Sex = Female      Weight = 155  
Test = DUI      Cit = 01100000000000  
Dr. Lic. = ND/NA  
Lot No = 01923080A3  
Cyl No = 23  
Expiration Date = 02/05/2025  
County = 09      Oper No. = 090801

  
\_\_\_\_\_  
Operator Signature  
KATHRYN ALLEN

Remarks:

Form 106-I8000

*Reinstad at Jail!*