



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59281 (06/2018)

| | |
|--|---|
| Serial Number 80-006492 | Instrument Location Corrections File room |
| Reason for Install/Repair | |
| <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____ | |

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

| | |
|-------------------------------|------------------------|
| Field Inspector Signature | Date 8/25/23 |
|-------------------------------|------------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | |
|---------------------------|------------------------------------|
| Reviewed/Certified By | Certified Date 29Aug2023 |
|---------------------------|------------------------------------|

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006492
Location = FYPD 8164.14.00 09/16
08/25/2023 15:10

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 14323080A1
Standard Cyl #? 17
Standard Expiration? 06/05/2025
Oper No? 130688

Flow Cal. Date: 08/18/2015
Slope 669
Intercept -631515

IR Calibration Date: 08/18/2015
 3um 9um

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -41685 | -23145 |
| 1st Coef(*100): | 269226 | 142194 |
| 2nd Coef(*100): | 2671 | 1239 |
| H2O adj(mg/l*10k): | 872 | 437 |

***** Printer Test End *****



Operator Signature
MONICA PATNEAUD

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006492
Location = FYPD 8164.14.00 09/16
08/25/2023 15:11

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:12 |
| 02 Std. Gas | 0.079 | 15:13 |
| 03 Room Air | 0.000 | 15:13 |
| 04 Std. Gas | 0.078 | 15:14 |
| 05 Room Air | 0.000 | 15:14 |
| 06 Std. Gas | 0.079 | 15:15 |
| 07 Room Air | 0.000 | 15:15 |


Lot No = 14323080A1

Cyl No = 17

Exp Date = 06/05/2025

County = 61

Oper No. = 130688



Operator Signature
MONICA PATNEAUD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006492
Location = FYPD 8164.14.00 09/16
08/25/2023 15:16

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:18 |
| 02 *Subject Test | RFI* | 15:18 |
| 03 Room Air | 0.000 | 15:18 |

*Invalid Test
Inhibited - RFI

Sub Name = JOHNSON, AMBER TEST
Sub DOB = 08/03/1984
Sub Sex = Female Weight = 150
Test = DUI Cit = TEST
Dr. Lic. = ND/NA
Lot No = 14323080A1
Cyl No = 17
Expiration Date = 06/05/2025
County = 61 Oper No. = 130688

*RFI
test
for
install*



Operator Signature
MONICA PATNEAUD

Remarks: