

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

Serial Number 80-000492 Instrument Location COV/Poth on 3 File VIDONA
Reason for Install/Repair
Install After Receiving From Crime Laboratory Install After Location Change
Other (Specify)
Check When Done:
1. Surge Protector Installed/Property Grounded.
2. Telephone Line Connected to Intoxilyzer® 8000.
3. Breath Tube Healed.
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
6/ Run Tests:
A. Print Test (Level1, Function P).
B. ACA Test (Level 1, Function C).
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
3. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.
Field Inspector Signature
0/45/23
Crime Laboratory Use Only
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.
O CONTRACTOR OF THE CONTRACTOR

Reviewed/Certified By

Certified Date 29Aug2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006492
Location = FYPD 8164.14.00 09/16
08/25/2023 15:10

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer (Y/N)? No No Display Volume? Yes Disable On Memfull? 1 # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 14323080A1 17 06/05/2025 Standard Expiration? 130688 Oper No?

Flow Cal. Date: 08/18/2015 Slope 669 Intercept -631515

IR Calibration Date: 08/18/2015 3um 9um

Oth	Coef(*100):	-41685	-23145
	Coef(*100):	269226	142194
2nd	Coef(*100):	2671	1239
H20	adj(mg/l*10k	:): 872	437

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature MONICA PATNEAUD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006492 Location = FYPD 8164.14.00 09/16 08/25/2023

15:11

## DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	15:12
02	Std.	Gas	0.079	15:13
03	Room	Air	0.000	15:13
04	Std.	Gas	0.078	15:14
05	Room	Air	0.000	15:14
06	Std.	Gas	0.079	15:15
07	Room	Air	0.000	15:15

Lot No = 14323080A1

Cyl No = 17

Exp Date = 06/05/2025

County = 61

Oper No. = 130688

Operator Signature MONICA PATNEAUD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FYPD
08/25/2023

SN 80-006492 8164.14.00 09/16 15:16

Test	AC	Time
01 Room Air	0.000	15:18
02 *Subject Test	RFI*	15:18
03 Room Air	0.000	15:18

\*Invalid Test Inhibited - RFI

Sub Name = JOHNSON, AMBER TEST

Sub DOB = 08/03/1984

Sub Sex = Female

Weight = 150 Cit = TEST

Test = DUI

Dr. Lic. = ND/NA Lot No = 14323080A1

Cyl No = 17

Expiration Date = 06/05/2025

County = 61

Oper No. = 130688

Operator Signature MONICA PATNEAUD

Remarks:

Form 106-I8000

RF1 Jest Jestall