

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 60 5957 Instrument Location Mandan					
Reason for Install/Repair					
Install After Receiving From Crime Laboratory Install After Location Change					
Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):	,				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	e Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the A					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair 104-G), Print Test, ACA Test, and RFI Test.					
Field Inspector Signature	Data				
The state of the s	9/6/0002				
Crime Laboratory Use Only	1 44 000				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to de	termina alcohal assesstation (
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	rime Laboratory Division, in the				
Reviewed/Certified By ULL	Certified Date 07Scp2023				

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005957
Location = MAND 8164.14.00 09/16
09/06/2023 15:57

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 14323080A1 73 06/05/2025 Standard Expiration? Oper No? 020611

Flow Cal. Date: 06/21/2013 Slope 647 Intercept -613317

IR Calibration Date: 05/11/2015
3um 9ur

Oth Coef(*100): -16666 -18045
1st Coef(*100): 278478 132949
2nd Coef(*100): 3303 1397
H2O adj(mg/l*10k): 379 351

******* Printer Test End ********

Operator Signature
TOM SHARP

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MAND
09/06/2023

Alcohol Analyzer SN 80-005957 8164.14.00 09/16 15:58

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:58
02	Std.	Gas	0.078	15:59
03	Room	Air	0.000	15:59
04	Std.	Gas	0.079	16:00
05	Room	Air	0.000	16:00
06	Std.	Gas	0.079	16:01
07	Room	Air	0.000	16:01

Lot No = 14323080A1

Cyl No = 73

Exp Date = 06/05/2025

County = 30

Oper No. = 020611

Operator Signature TOM SHARP

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = MAND 09/06/2023

SN 80-005957 8164.14.00 09/16 16:02

DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	16:03
02	Room	Air	RFI*	16:03

*Invalid Test Inhibited - RFI

Lot No = 14323080A1

Cyl No = 73

Exp Date = 06/05/2025

County = 30

Oper No. = 020611

Operator Signature TOM SHARP

Remarks:

Form 106-I8000