

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Control Number 173				
Serial Number 80-003061 Instrument Location 6FCCC				
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
✓/3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair 104-G), Print Test, ACA Test, and RFI Test.	Checkout (SFN59281, Form			
Field Impector Signature	Date 9/8/2023			
rime Laborator Use Only	7/			
This installation has been regiewed and the instrument is approved to be used for the analysis of breath to det the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crounty of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	ermine alcohol concentration from rime Laboratory Division, in the			
Reviewed/Certified By Mul	Certified Date 08Sep2023			

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = GFSO 8164.14.00 09/16
09/08/2023 07:43

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? NoDisplay Volume? $N \circ$ Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 18 02/05/2025 Standard Expiration? 020836 Oper No?

Flow Cal. Date: 03/04/2011 Slope 700 Intercept -623726

IR Calibration Date: 04/10/2015 3um 9um

Oth Coef(*100): -25722 -17879

1st Coef(*100): 266163 136703

2nd Coef(*100): 2302 1151

H2O adj(mg/l*10k): 412 294

****** Printer Test End *******

Oper t/ Signatu/ D/ D STROMBE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFSO
09/08/2023

Alcohol Analyzer SN 80-003061 8164.14.00 09/16 07:46

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	07:47
02	Std.	Gas	0.079	07:47
03	Room	Air	0.000	07:48
04	Std.	Gas	0.080	07:48
05	Room	Air	0.000	07:49
06	Std.	Gas	0.080	07:49
07	Room	Air	0.000	07:50

Lot No = 01923080A3

Cyl No = 18

Exp Date = 02/05/2025

County = 18

Oper No. = 020836

erator Signatu DAVID STROMBER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = GFSO 09/08/2023

SN 80-003061 8164.14.00 09/16 07:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	07:45
02 Reference	RFI*	07:45
03 Room Air	0.000	07:45

*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 18

Exp Date = 02/05/2025

County = 18

Oper No. = 020836

ator Sign ure AVID STROWERG

Remarks:

Form 106-I8000