

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number			
80-006508	Instrument Location		
Reason for Install/Repair	RICHLAND CO. JAZL; WAHPETON		
Install After Receiving From Crime Laboratory Install After Location Change			
& Other (Specify) Annual Inspection a	+ Chine Laboratory_		
Check When Done:	0		
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	tion; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
🔀 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
	Alc		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form		
Field This pector Signature			
Chilles all	Date		
Anyton 10m	06/06/2013		
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to b	e used for the analysis of breath to determine allochal approximation (

the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	
hand	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006508 8164.14.00 09/16 Location = WAHP 17:00 06/06/2023 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 34921080A1 35 02/05/2024 Standard Expiration? 132570 Oper No? Flow Cal. Date: 08/20/2015 644 Slope -568181 Intercept IR Calibration Date: 08/20/2015 3um 9um ______ -18627 0th Coef(*100): -22365 1st Coef(*100): 256969 2nd Coef(*100): 3168 136095 1444 453 H2O adj(mg/l*10k): 693 ********** Printer Test End **********

Operator Signature CHRISTOPHER HUARD


Remarks: PRENT TOST

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006508Location = WAHP8164.14.00 09/1606/06/202316:56

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.079 0.000 0.078 0.000	16:56 16:57 16:57 16:58 16:58 16:58 16:59



Derator Signature CHRISTOPHER HUARD

Remarks: ACA TEST

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer SN 80-006508 North Dakota Model 8000 8164.14.00 09/16 Location = WAHP 17:01 06/06/2023 Test AC Time 17:02 01 Room Air RFI* 17:03 02 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = DOE, JOHN NA Sub DOB = 01/01/1980Sub Sex = Male Weight = 215Cit = 123456789Test = DUI Dr. Lic. = ND/DOE803432Lot No = 34921080A1Cyl No = 35Expiration Date = 02/05/2024Oper No. = 132570County = 39

Operator Signature CHRISTOPHER HUARD

Remarks: **RF***I* TEST

Form 106-18000