



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-006687</b>	Instrument Location <b>PLPD</b>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Annual Inspection at Crime Laboratory</u>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

*uw  
10JUL23*

Field Inspector Signature <i>[Signature]</i>	Date <b>07-03-2023</b>
---	---------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <b>10 July 2023</b>
---	---------------------------------------

*No tests ran after last chain of custody, so no install paperwork was needed. AEN 10Jul23*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006687  
Location = PLPD                    8164.14.00 09/16  
06/30/2023                            11:21

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34921080A1  
Standard Cyl #?            34  
Standard Expiration?        02/05/2024  
Oper No?                    132720

Flow Cal. Date:            10/05/2018  
Slope                        696  
Intercept                   -609977

IR Calibration Date:        10/05/2018

	3um	9um
0th Coef(*100):	-11478	-18771
1st Coef(*100):	255256	132135
2nd Coef(*100):	4432	1676
H2O adj(mg/l*10k):	604	575

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
BRYTON ALLARD

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006687  
Location = PLPD      8164.14.00 09/16  
06/30/2023      11:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:25
02 Std. Gas	0.079	11:25
03 Room Air	0.000	11:26
04 Std. Gas	0.079	11:26
05 Room Air	0.000	11:27
06 Std. Gas	0.078	11:27
07 Room Air	0.000	11:27

Lot No = 34921080A1  
Cyl No = 34  
Exp Date = 02/05/2024  
County = 07      Oper No. = 132720



Operator Signature  
BRYTON ALLARD

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006687  
Location = PLPD                8164.14.00 09/16  
06/30/2023                      11:29

Test	AC	Time
01 Room Air	RFI*	11:30
02 Room Air	RFI*	11:30

\*Invalid Test  
  Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 01/01/2000  
Sub Sex = Male                    Weight = NA  
Test = OTH                        Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 34921080A1  
Cyl No = 34  
Expiration Date = 02/05/2024  
County = 07                        Oper No. = 132720



Operator Signature  
BRYTON ALLARD

Remarks:

Form 106-I8000