

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	nstrument Location	day this old	
Reason for Install/Repair	SATGENT county 51	HAUM'S OFFICE.	
Install After Receiving From Crime Laboratory Install After L	Location Change		
Other (Specify)	. •		
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.	a.		
3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Lev	vel 1, Function B or C; Key Radio Du	ring Test).	
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5049	6, Form 120-G) and Place it by the I	ntoxilyzer® for Use.	
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Intoxily 104-G), Print Test, ACA Test, and RFI Test.	yzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature		Date	
Ja word		11/03/23	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be us the date the Field Inspector performed the installation. This record on file at County of Burleigh, North Dakota, is certified to be a true and correct copy of	t the Office of Attorney General, Crin		
Reviewed/Certified By		Certified Date 15Nov2023	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006690
Location = SACO 8164.14.00 09/16
11/03/2023 11:40

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34921080A1 Standard Cyl #? 57 Standard Expiration? 02/05/2024 Oper No? 132607

Flow Cal. Date: 08/23/2018 Slope 677 Intercept -593298

IR Calibration Date: 08/23/2018
3um 9um

3um 9um

Oth Coef(\*100): -8380 -21663

1st Coef(\*100): 265142 130429

2nd Coef(\*100): 3933 1429

H2O adj(mg/l\*10k): 412 539

\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
JONATHAN WOLD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = SACO
11/03/2023

Alcohol Analyzer SN 80-006690 8164.14.00 09/16 11:34

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:36
02 Std. Gas	0.081	11:36
03 Room Air	0.000	11:37
04 Std. Gas	0.081	11:37
05 Room Air	0.000	11:38
06 Std. Gas	0.081	11:38
07 Room Air	0.000	11:39

Lot No = 34921080A1

Cyl No = 57

Exp Date = 02/05/2024

County = 41

Oper No. = 132607

Operator Signature
JONATHAN WOLD

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = SACO
11/03/2023

Alcohol Analyzer SN 80-006690 8164.14.00 09/16 11:43

## DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	11:44
02 Room Air	0.000	11:44

\*Invalid Test Inhibited - RFI

Lot No = 34921080A1

Cyl No = 57

Exp Date = 02/05/2024

County = 41

Oper No. = 132607

Operator Signature JONATHAN WOLD

Remarks:

Form 106-I8000