

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location	
Reason for Install/Repair	
Install After Receiving From Crime Laboratory Install After Location Change	
Prother (Specify) vpdute	
Check When Done:	
Z 1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
☐ 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Ir	
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form
Field Inspector Signature	Date 12/15/23
Crime Laboratory Use Only	/ -/ 2
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to deter the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crim County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	mine alcohol concentration from ne Laboratory Division, in the

Reviewed/Certified By Au All Certified Date 14Dec2023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004944 Location = WILL 8164.14.00 09/16 12/05/2023 15:11 abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A4 Standard Cyl #? 43 Standard Expiration? 06/05/2025 Oper No? 130739 Flow Cal. Date: 07/01/2011 Slope 685 Intercept -634740 IR Calibration Date: 07/01/2011 3um 9um _____ 0th Coef(*100): -22591 -172341st Coef(*100): 265915
2nd Coef(*100): 2841 132691 1562 H2O adj(mg/l*10k): 559425 ********* Printer Test End *********

Operator Signature MICHELLE ROMANS

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004944Location = WILL8164.14.00 09/1612/05/202315:12

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.083 0.000 0.082 0.000 0.082 0.000	15:12 15:13 15:13 15:13 15:14 15:14 15:15

Lot No = 14323080A4 Cyl No = 43 Exp Date = 06/05/2025 County = 53 Oper No

Oper No. = 130739

Operator signature MICHELLE ROMANS

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-004944 Location = WILL 8164.14.00 09/16 12/05/2023 15:15 Test AC Time 01 Room Air RFI* 15:17 02 Room Air 0.000 15:17 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 01/19/1990Sub Sex = Male Weight = 123Test = DUI Cit = TEST Dr. Lic. = ND/TEST Lot No = 14323080A4Cyl No = 43Expiration Date = 06/05/2025 County = 53Oper No. = 130739

Operator Signature MICHELLE ROMANS

Remarks:

Form 106-18000

. •