

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-00 6000 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No.

- B. Calibration Adjustment (Level 3, M, C, O)
 - 1. Autocalibration Printout Attached
 - Max Power Res Value ≥ 10
 - Auto Range Res Value ≥ 4
 - 2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H ₂ O	NA-Milli-Q H ₂ O	MP5317
2	0.040	202303H	28 Maras	MP5320
3	0.080	202302B	14 Febas	MP5290
4	0.100	202304A	OHApras	MP3068
5	0.300	202201F	19Jan 24	MP3066

- 3. 0.080 AC Calibration Gas for H₂O Adjustment Lot No. <u>IH323080A4</u> Cyl No. <u>013</u> Exp. Date: <u>6|5|25</u>
- 4. Atmospheric Pressure
 - Displayed by Intoxilyzer[®] 8000 Adjusted to using barometer Auto Calibration Report printout Barometer Model Barometer Serial Number Barometer Calibration Expiration Date

<u>956</u> mbar <u>966</u> mbar <u>965</u> mbar <u>10510-922</u> <u>230307250</u> 02 May2025

- 5. Screen displayed "Calibration Success"
- 6. X Calibration Adjustment Printout Attached
 ➢ Solution 1 Avg % Abs ≤ 0.2500
 Solution 2-5 REL STD DEV ≤ 3.000

Intoxilyzer 8000 Calibration Adjustment Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED Qualtrax ID: 11859 Revision: 1 Status: Published Date Approved: 06/20/2023 Page 1 of 2 Residual (g/210 L) values for solutions 1 - 5 \leq 0.0020 for 3 µm and 9 um channels Dry Gas H₂O adjustment sum for 3 µm and 9 µm channels within ± 10 3 μm <u>3511</u> (Ave.) + <u>232</u> (H₂O Adj.) = <u>3809</u> 9 µm <u>3432</u> (Ave.) + <u>31</u> (H₂O Adj.) = <u>3809</u>

C. Is an Annual Inspection due for this instrument? Yes or No If Yes, complete Intoxilyzer 8000 Annual Inspection (Qualtrax ID: 11698) If No, complete Intoxilyzer 8000 Calibration (Qualtrax ID: 11871).

Remarks/Notes:

alyst Signature

C Pusschaller eviewer Signature

12 Dec 2023

Intoxilyzer 8000 Calibration Adjustment Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

Qualtrax ID: 11859 Revision: 1 Status: Published Date Approved: 06/20/2023 Page 2 of 2 TOXL Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006666 12/12/2023 13:56:02

Auto Calibration

<<<< 3um >>>> <<<< 9um >>>>

 Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1

 Sample
 % Abs
 (% Abs Ref)
 % Abs
 (% Abs Ref)

 Sample #1
 0.0660
 (0.0040)
 0.2430
 (-0.0020)

 Sample #2
 0.0940
 (0.0240)
 0.2460
 (0.0000)

 Sample #3
 0.0860
 (0.0610)
 0.2330
 (0.0130)

 Sample #4
 0.0650
 (0.0840)
 0.2360
 (0.0210)

 Avg % Abs
 0.0817
 (0.0563)
 0.2383
 (0.0113)

 STD DEV
 0.0150
 (0.0303)
 0.0068
 (0.0106)

 REL STD DEV
 18.340
 (53.736)
 2.856
 (93.518)

 Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1

 Solution = 0.040 g/210L or 0.1905 mg/1, Samples = 4, Discarded = 1

 Sample
 % Abs
 (% Abs Ref)
 % Abs
 (% Abs Ref)

 Sample #1
 0.7930
 (-0.0200)
 1.6360
 (-0.0100)

 Sample #2
 0.7980
 (-0.0040)
 1.6170
 (0.0120)

 Sample #3
 0.7630
 (0.0230)
 1.5870
 (0.0350)

 Sample #4
 0.7710
 (0.0230)
 1.5890
 (0.0330)

 Avg % Abs
 0.7773
 (0.0140)
 1.5977
 (0.0267)

 STD DEV
 0.0183
 (0.0156)
 0.0168
 (0.0127)

 REL STD DEV
 2.359
 (111.346)
 1.050
 (47.779)

 Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1

 Solution = 0.080 g/210L or 0.3810 mg/1, Samples = 4, Discarded = 1

 Sample
 % Abs
 (% Abs Ref)
 % Abs
 (% Abs Ref)

 Sample #1
 1.5190
 (-0.0210)
 2.9410
 (-0.0100)

 Sample #2
 1.5260
 (0.0010)
 2.9640
 (-0.0100)

 Sample #3
 1.4890
 (0.0220)
 2.9480
 (-0.0100)

 Sample #4
 1.4980
 (0.0290)
 2.9590
 (-0.0120)

 Avg % Abs
 1.5043
 (0.0173)
 2.9570
 (-0.0077)

 STD DEV
 0.0193
 (0.0146)
 0.0082
 (0.0059)

 REL STD DEV
 1.283
 (84.067)
 0.277
 (76.428)

 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

 Sample
 % Abs
 (% Abs Ref)
 % Abs
 (% Abs Ref)

 Sample #1
 1.8570
 (-0.0050)
 3.5930
 (-0.0020)

 Sample #2
 1.8530
 (0.0220)
 3.6050
 (0.0150)

 Sample #3
 1.8430
 (0.0120)
 3.6140
 (0.0060)

 Sample #4
 1.8310
 (0.0260)
 3.6050
 (0.0220)

 Avg % Abs
 1.8423
 (0.0200)
 3.6080
 (0.0143)

 STD DEV
 0.0110
 (0.0072)
 0.0052
 (0.0080)

 REL STD DEV
 0.598
 (36.056)
 0.144
 (55.959)

 Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

 Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

 Sample
 % Abs
 (% Abs Ref)
 % Abs
 (% Abs Ref)

 Sample #1
 5.0830
 (-0.0210)
 9.6080
 (-0.0040)

 Sample #2
 5.1470
 (-0.0020)
 9.7090
 (0.0250)

 Sample #3
 5.1620
 (0.0070)
 9.7350
 (0.0420)

 Sample #4
 5.1820
 (-0.0050)
 9.7050
 (0.0520)

 Avg % Abs
 5.1637
 (0.0000)
 9.7163
 (0.0397)

 STD DEV
 0.0176
 (0.0062)
 0.0163
 (0.0137)

 REL STD DEV
 0.340
 (8046618112.000)
 0.168
 (34.413)

pg 1 of 2

TOXL Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006666 12/12/2023 13:56:02

Auto Calibration

< •	<<<< 3ur	n >>>>	<<<<<	9um :	>>>>>
Zero Order Coef -197.99 First Order Coef 2628.72 Second Order Coef 34.08				4.09 1.57 3	
Act (g/210L) 0.000 0.040 0.080 0.100 0.300	Fit (g/210L) 0.000 0.039 0.081 0.100 0.300	Residual (g/210L) -0.0004 0.0008 -0.0005 0.0000 0.0000	Act (g/210L) 0.000 0.040 0.080 0.100 0.300	Fit (g/210L 0.000 0.040 0.080 0.100 0.300	Residual) (g/210L) -0.0002 0.0004 -0.0002 -0.0001 0.0000

	<<<<<	3um	>>>>>	<<<<<	9um >>>>
Solution = Sample	0.080 g/21	0L or 0	.3810 mg/l,	Samples = 4,	Discarded = 1
Sample #1		35'	71.00		3418.00
Sample #2		350	53.00		3435.00
Sample #3		350	59.00		3429.00
Sample #4		36	00.00		3432.00
Avq		35'	77.3333		3432.0000
STD DEV		19	.8578		3.0000
REL STD DEV		0.	555		0.087
H2O adjust	(mg/l*10k)	23:	2		377

Atmospheric Pressure = 965

TOXL				
intoxilyzer	-	Alcohol	Analyzer	
Model 8000			SN	80-005865
12/12/2023				13:56:92

Auto Calibration Max Power Res Value = 47 Auto Range Res Value = 24 pg 2 of 2