

INTOXILYZER® 8000 CALIBRATION

Intoxilyzer [®] 8000 Serial Number: 80-00 6444 Calibration Location: TOXL							
Α.	A. Pre-Calibration						
	1.	🛛 Download Data					
	2.	KUpload Operator File					
	3.	Kattery Check	\sim				
		Was the external battery pack replaced?	Yes or No				
	4.	O-Rings	Ver er D				
		Replaced Simulator O-Ring Replaced Breath Tube O-Ring	Yes or No Yes or No				
		Replaced Breath Tube O-Ring	les of No				
В.	Gener	al Şetup and Checks:					
	1.	Diagnostics passed and instrument in "Ready" mod	le				
	2.	Breath tube heated					
	3.	Date, time and location code (Level 2,E). Re-set if					
		Time Zone: CST or MDT (Time on test records	s will be in time zone circled).				
	4.	Print test (Level 1,P). Sign and attach test record.					
	5.	Tank monitor (Level 3,D,G).	nei				
	Display: <u>913</u> psi Regulator: <u>915</u> psi Display and Regulator ± 50 psi of each other ves or No						
	Gas tank tare necessary? Yes or No						
		If Yes, display readings after tare (Level 3,M,C					
		Display: psi Regulator:	psi				
C	Q. Tasta (Qinn and attach tast records):						
Ο.	1.	(Sign and attach test records): ☑ Configure simulator for the following test (Level 1,5)	3)				
	1.	Wet Calibration Check - Low AC (Level 1,C)	5).				
	Known Value $\leq 0.03 \text{ AC}$: 0.03 D AC						
		Sim. Ser #: <u>MP3D65</u>					
	Lot #:						
	Exp. Date: <u>IJana4</u>						
		Results ± 0.005 of known AC					

Intoxilyzer 8000 Calibration

Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director

UNCONTROLLED WHEN PRINTED

Qualtrax ID: 11871 Revision: 2 Status: Published Date Approved: 06/21/2023 Page 1 of 2 2. \bigtriangledown Configure simulator for the following test (Level 1,S). Wet Calibration Check - High AC (Level 1,C) Known Value $\geq 0.25 \text{ AC}$: <u>0.400</u> AC Sim. Ser #: MP3059 Lot #: 2022/2C Exp. Date: 08Dec24 Results ± 5% AC of known AC Configure dry gas standard for the remaining tests (Level 1,S). 3. Known Value : 0.080 AC Gas Cylinder Lot #: 14323080A 4 Cylinder #: 013 Exp. Date: 6525 4. X Interferent Check (Level 1,B) Known Value: 0.10 AC + 0.05% Acetone Lot #: ICST Exp. Date: NIA Display reads "Interferent Detect" XRFI Check (CMS Mode) 5. Display reads "RFI Detect" Dry Calibration Check (Level 1,C) 6. Test 1 0.079 Test 4 0.079 Test 7 0,079 Test 8 0.079 Test 5 0.018 Test 2 0.079 Test 9 0.019 Test 3 0.079 Test 6 0.019 Average 0.018

Results ± 0.005 AC of known AC

D. Remarks/Maintenance: NHA Completed after a Calibration justment was completed

Instrument is acceptable to be used in the field. (Yes or No

If No, state reason(s) why:

alibration Analyst Signature

dle Portschuller viewer

Intoxilyzer 8000 Calibration Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

<u>15Dec2023</u> Date <u>15 Dec 2023</u>

Qualtrax ID: 11871 Revision: 2 Status: Published Date Approved: 06/21/2023 Page 2 of 2

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006666 Location = TOXL 8164.14.00 09/16 12/15/2023 12:47 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A4 Standard Cyl #? 13 Standard Expiration? 06/05/2025 Oper No? 133237 Flow Cal. Date: 05/16/2023 Slope 686 Intercept -761238 IR Calibration Date: 12/12/2023 3um 9um -------Oth Coef(*100): -19799 1st Coef(*100): 262872 2nd Coef(*100): 3408 -31409 135157 1553 H2O adj(mg/l*10k): 232 377 ********** Printer Test End **********

Operator Signature ANNA NAREHOOD

Remarks:

Print Test

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006666Location = TOXL8164.14.00 09/1612/15/202313:04

		WET CAL CHECK	
Τe	est	AC	Time
01	Room Air	0.000	13:04
02	Std. Sol	. 0.028	13:05
03	Room Air	0.000	13:06
04	Std. Sol	. 0.028	13:06
05	Room Air	0.000	13:07
06	Std. Sol	. 0.028	13:08
07	Room Air	0.000	13:08

08 Sim Temp = 34.0°C

Simul Ser No = MP3065 Std Sol No = 202201B County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

OW AC Check-0.030AC

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006666 8164.14.00 09/16 Location = TOXL 12/15/2023 13:22

WET CAL CHECK AC 01 Room Air 0.000

Time

02 03 04 05	Room Std. Room Std. Room Std.	Sol. Air Sol. Air		13:23 13:25 13:25 13:26 13:27 13:27
06	Std.	Sol.	0.406	13:27
07	Room	Air	0.000	13:28

08 Sim Temp = 34.0°C

Simul Ser No = MP3059 Std Sol No = 202212COper No. = 133237 County = 08

Operator Signature ANNA NAREHOOD

Remarks:

Test

High AC Check-0.400AC

NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006666 Location = TOXL 8164.14.00 09/16 12/15/2023 13:30 Test Time AC 01 Room Air 0.000 13:31 02 *Subject Test INT* 13:31 03 Room Air 0.000 13:32 *Invalid Test Interferent Detected Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Test = OTHCit = INTERFERENT CK Dr. Lic. = ND/TES989643 Lot No = 14323080A4Cyl No = 13Expiration Date = 06/05/2025County = 08Oper No. = 133237

Intoxilyzer Test Record and Checklist

Operator Signature ANNA NAREHOOD

Remarks: Interferent Cheuk

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006666 8164.14.00 09/16 Location = TOXL 12/15/2023 13:32 AC Time Test 13:33 01 Diagnostic OK 13:34 RFI* 02 Room Air 13:34 03 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Cit = RFI CHECK Test = OTHDr. Lic. = ND/TES989643Lot No = 14323080A4Cyl No = 13Expiration Date = 06/05/2025County = 08Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks:

RFI Check

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006666Location = TOXL8164.14.00 09/1612/15/202313:34

DRY CAL CHECK Time AC Test 01 Room Air 0.000 13:35 02 Std. Gas 0.079 13:35 03 Room Air 0.000 13:36 04 Std. Gas 0.079 13:36 13:36 05 Room Air 0.000 06 Std. Gas 0.079 13:37 07 Room Air 0.000 13:37

Lot No = 14323080A4 Cyl No = 13 Exp Date = 06/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Cal cherk #1-3

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006666Location = TOXL8164.14.00 09/1612/15/202313:38

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	13:38
02	Std.	Gas	0.079	13:39
03	Room	Air	0.000	13:39
04	Std.	Gas	0.078	13:39
05	Room	Air	0.000	13:40
06	Std.	Gas	0.079	13:40
07	Room	Air	0.000	13:41

Lot No = 14323080A4 Cyl No = 13 Exp Date = 06/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Chey #4-6

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006666
Location = TOXL	8164.14.00 09/16
12/15/2023	13:41

			DRY CAL CHECK	
Т	est		AC	Time
01	Room	Air	0.000	13:42
02	Std.	Gas	0.079	13:42
03	Room	Air	0.000	13:42
04	Std.	Gas	0.079	13:43
05	Room	Air	0.000	13:43
06	Std.	Gas	0.079	13:44
07	Room	Air	0.000	13:44

Lot No = 14323080A4 Cyl No = 13 Exp Date = 06/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Dry cal check #7-9