



ETHANOL GAS STANDARD CYLINDER REPORT
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (08/2018)

| | | |
|--|---|---------------------------------------|
| Chemical Test Operator Name (Print) <i>Christopher N. Ormiston</i> | | |
| Location <i>CARR</i> | Intoxilyzer® Serial Number <i>80-006504</i> | |
| Gas Lot Number <i>01923080A3</i> | Gas Cylinder Number <i>25</i> | Gas Expiration Date <i>02/05/2025</i> |

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. *080* AC
 - B. 0. *081* AC
 - C. 0. *080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| | |
|--|------------------------|
| Chemical Test Operator Signature <i>[Signature]</i> 1621 | Date <i>12/15/2023</i> |
| Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i> | Date <i>29Dec2023</i> |


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006504
Location = CARR 8164.14.00 09/16
12/15/2023 13:50

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 13:51 |
| 02 Std. Gas | 0.080 | 13:51 |
| 03 Room Air | 0.000 | 13:52 |
| 04 Std. Gas | 0.081 | 13:52 |
| 05 Room Air | 0.000 | 13:52 |
| 06 Std. Gas | 0.080 | 13:53 |
| 07 Room Air | 0.000 | 13:53 |

Lot No = 01923080A3
Cyl No = 25
Exp Date = 02/05/2025
County = 16 Oper No. = 131530


Operator Signature
CHRISTOPHER ORMISTON

Remarks: ACA Test