

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	0			
Serial Number 80 - 004950	Instrument Location BMDC	יני ט א		
Reason for Install/Repair	Barrenge County ster	to Dept.		
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change			
Other (Specify)	er Education Griange			
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	ion; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
☐ 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).		
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date 1/17/24		
Crime Laboratory Use Only		ι /		
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from				
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By		Certified Date		
		Certified Date		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = BMDC 8164.14.00 09/16
01/17/2024 09:57

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 01923080A3 53 Standard Expiration? 02/05/2025 Oper No? 131854

Flow Cal. Date: 01/18/2010 Slope 774 Intercept -908904

IR Calibration Date: 05/05/2016

3um 9um

Oth Coef(\*100): -13717 -19495
1st Coef(\*100): 271916 140513
2nd Coef(\*100): 3504 1503
H2O adj(mg/l\*10k): 370 426

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature NICK SCHOEDEL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = BMDC 01/17/2024

SN 80-004950 8164.14.00 09/16 10:04

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:05
02 Std. Gas	0.084	10:06
03 Room Air	0.000	10:06
04 Std. Gas	0.083	10:06
05 Room Air	0.000	10:07
06 Std. Gas	0.082	10:07
07 Room Air	0.000	10:08

Lot No = 01923080A3

Cyl No = 53

Exp Date = 02/05/2025

County = 08

Oper No. = 131854

Operator Signature NICK SCHOEDEL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Al North Dakota Model 8000 Location = BMDC 81 01/17/2024

Alcohol Analyzer SN 80-004950 8164.14.00 09/16 10:12

Τe	est	AC	Time
	Room Room	 RFI* 0.000	10:14

\*Invalid Test Inhibited - RFI

Sub Name = TEST, RFI NA Sub DOB = 01/01/2001

Sub Sex = Male Weight = NA Test = OTH Cit = NA

Dr. Lic. = ND/NA
Lot No = 01923080A3
Cyl No = 53

Expiration Date = 02/05/2025

County = 08 Oper No. = 131854

Operator Signature
NICK SCHOEDEL

Remarks:

RFI Yest no Subject

Form 106-I8000