

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

80-006678	Instrument Location GFP0			
Reason for Install/Repair				
🔀 Install After Receiving From Crime Laboratory 🛛 Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
- 🛱 3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
🖾 6. Run Tests:				
A. Print Test (Level1, Function P).				
🕅 B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Dur	ing Test).		
7. Repair and/or Maintenance Performed (if any): uplatez	user list at crime !	ab		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
步9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)	) at the Intoxilyzer® Location at the Ager	ıcy.		
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Ch	neckout (SFN59281, Form		
Field Inspector Signature		Date		
Aton Xann		OF-2 Fay		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	0	2	0	
Reviewed/Certified By	Am	Mal	$\rho$	Certified Date 29Jan2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006678 Location = GFPD8164.14.00 09/16 01/27/2024 14:30 abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 57 Standard Expiration? 02/05/2025 039947 Oper No? 08/22/2016 Flow Cal. Date: Slope 686 Intercept -611596 08/22/2016 IR Calibration Date: 3um 9um \_\_\_\_\_\_ \_\_\_\_\_ 0th Coef(\*100): -15445 -21650 1st Coef(\*100): 251854 131766 2nd Coef(\*100): 3094 1378 H2O adj(mg/l\*10k): 549513 \*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*

039947

Operator Signature JONATHAN LAMPI

Remarks:

Install after repair at crime lab

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006678Location = GFPD8164.14.00 09/1601/27/202414:30

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.083 0.000 0.083 0.000 0.083 0.000	14:31 14:32 14:32 14:33 14:33 14:33 14:33

039947 and

Operator Signature JONATHAN LAMPI

Remarks:

Install after repair at crime lab,

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006678Location = GFPD8164.14.00 09/1601/27/202414:35

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Room Air	0.000 0.083 RFI* 0.000	14:35 14:36 14:36 14:36
*Invalid Test Inhibited -	RFI	

39947

Operator Signature JONATHAN LAMPI

Remarks:

Install afterrepair at crimelab

Form 106-18000

,