

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-004944	Williston, ND	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory	Install After Location Change	
Other (Specify)	-	

Check When Done:

1 Surga Protector Installed/Property Grounded	
✓ 1. Surge Protector Installed/Property Grounded.	
$\square$ 2. Telephone Line Connected to Intoxilyzer® 8000.	
Z 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio D	uring Test).
7. Repair and/or Maintenance Performed (if any): N/A	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Ag	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair 104-G), Print Test, ACA Test, and RFI Test.	Checkout (SFN59281, Form
Field Inspector Signature	Date 22/07/2024
- V	

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	$\cap$	0		7	
Reviewed/Certified By	Ň	ne Mal	X		Certified Date 08Feb2024

and the

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004944 Location = WILL 8164.14.00 09/16 02/07/2024 17:45 abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABA Display Prelim Rslt? Yes DABACABA Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 # OI Princ Copies.Select Std (D/W/I)?DryStandard Value?0.080 Standard Lot #? Standard Cyl #? 14323080A4 43 Standard Expiration? 06/05/2025 Oper No? 132166 Flow Cal. Date: 07/01/2011 Slope 685 -634740 Intercept IR Calibration Date: 07/01/2011 9um 3um \_\_\_\_\_ 0th Coef(\*100): -22591 -17234 1st Coef(\*100): 265915 2nd Coef(\*100): 2841 132691 1562 H2O adj(mq/l\*10k): 559425 

Signa Dperator ALYSON MOSES

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004944
Location = WILL	8164.14.00 09/16
02/07/2024	17:46

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	17:46 17:47 17:47 17:47 17:48 17:48 17:48

Lot No = 14323080A4 Cyl No = 43 Exp Date = 06/05/2025 County = 53

Oper No. = 132166

e. Operat ALYSON MOSES

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004944 Location = WILL 8164.14.00 09/16 02/07/2024 17:50 Test Time AC 01 Room Air RFI\* 17:52 02 Room Air 0.000 17:52 \*Invalid Test Inhibited - RFI Sub Name = RFI TEST, INSTALL NA Sub DOB = 02/07/2004Sub Sex = Male Weight = NA Test = DUI Cit = NADr. Lic. = ND/NA Lot No = 14323080A4Cyl No = 43Expiration Date = 06/05/2025County = 53Oper No. = 132166

ALYSON MOSES

Remarks:

Form 106-18000