

#### NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

## **INTOXILYZER® 8000 ANNUAL INSPECTION**

Intoxilyzer	® 8000 Serial Number: 80-00 4948 Inspection Location: 10XL
A. Pre 1. 2. 3. 4. 5. 6.	Items with Instrument:   Gas Cylinder
B. Ge 1. 2. 3. 4. 5.	neral Setup and Checks:  ☐ Diagnostics passed and instrument in "Ready" mode ☐ Breath tube heated ☐ Date, time and location code (Level 2,E). Re-set if necessary.  ☐ Time Zone: CST or MDT (Time on test records will be in time zone circled). ☐ Print test (Level 1,P). Sign and attach test record. ☐ Tank monitor (Level 3,D,G). ☐ Display:
C. Te. 1.	sts (Sign and attach test records):  Configure simulator for the following test (Level 1,S).  Wet Calibration Check - Low AC (Level 1,C)  Known Value ≤ 0.03 AC: □ · ○ 2 · ○ AC  Sim. Ser #: NP3002  Lot #: 202301E  Exp. Date: □ 10005 of known AC

Intoxilyzer 8000 Annual Inspection

Document ID: 11698 Revision: 2

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Status: Published

**AEN** 

Approved By: Laboratory Director

Date Approved: 02/29/2024

UNCONTROLLED WHEN PRINTED Uploaded 05April2024

Page 1 of 2 Page 1 of 11

2.	Configure simulator for the following test (Level 1, Wet Calibration Check - High AC (Level 1, C)  Known Value ≥ 0.25 AC: 0.400 AC  Sim. Ser #: MP 8036  Lot #: 202212C  Exp. Date: 080ec24  X Results ± 5% AC of known AC	S).		
3.	Configure dry gas standard for the remaining tests	(Level 1,S).		
	Known Value: 0.080 AC Gas Cylinder Lot #: 14323080A4 Cylinder #: 13 Exp. Date: 0/5/25			
4.				
	Known Value: O.IOAC + 0.05% Acetone Sim. Ser #: DR7344			
	Lot #: 1CS 8			
	Exp. Date: NIA			
	☑ Display reads "Interferent Detect"			
5.	☐ RFI Check (CMS Mode)			
	☑ Display reads "RFI Detect"			
6.	Dry Calibration Check (Level 1,C)			
	Test 1 0.079 Test 4 0.079 Test 7 Test 2 0.078 Test 5 0.079 Test 8 Test 3 0.079 Test 6 0.079 Test 8 Average 0.078 Results ± 0.005 AC of known AC			
D. Rema	irks/Maintenance: Initially ran #1-3 dry cal c	heck (.084, .083,		
	so performed calibration adjustment			
	are all post-calibration adjustment			
	s acceptable to be used in the field. Yes or No			
f No, state reason(s) why:				
f Yes, chang	ge location code back to A.4. 💢			
		010002001		
Inspector Sig	gnature	OIApr2024 Date		
Reviewer	Portschiller	Od Apr 2024 Date		

Intoxilyzer 8000 Annual Inspection

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Approved By: Laboratory Director

UNCONTROLLED WHEN PRINTED Uploaded 05April2024

Document ID: 11698 Revision: 2

Status: Published Date Approved: 02/29/2024

Page 2 of 2

**AEN** 

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004948
Location = TOXL 8164.14.00 09/16
04/01/2024 13:17

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 21 Standard Expiration? 02/05/2025 133237 Oper No?

Flow Cal. Date: 08/28/2018

Slope 642 Intercept -606282

IR Calibration Date: 08/28/2018
3um 9um

Oth Coef(\*100): -13609 -14548 1st Coef(\*100): 248186 128726 2nd Coef(\*100): 3538 1508

H2O adj(mg/l\*10k): 602 515

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature ANNA NAREHOOD

Remarks: Print Test

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/01/2024

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 13:18

# WET CAL CHECK

Test		AC		Time
01	Room	Air	0.000	13:19
02	Std.	Sol.	0.021	13:19
03	Room	Air	0.000	13:20
04	Std.	Sol.	0.022	13:21
05	Room	Air	0.000	13:21
06	Std.	Sol.	0.022	13:22
07	Room	Air	0.000	13:22

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$ 

Simul Ser No = MP3002 Std Sol No = 202301E

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: LOW AC Check - 0.020AC

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = TOXL 04/01/2024

SN 80-004948 8164.14.00 09/16 13:24

#### WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:24
02 Std. Sol.	0.396	13:25
03 Room Air	0.000	13:26
04 Std. Sol.	0.396	13:26
05 Room Air	0.000	13:27
06 Std. Sol.	0.397	13:28
07 Room Air	0.000	13:28

08 Sim Temp = 34.0°C

Simul Ser No = MP6036 Std Sol No = 202212C

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

High Ac Check- 0.400AC Remarks:

SN 80-004948

13:44

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = TOXL 8164.14.00 09/16 04/01/2024

Test	AC	Time
01 Room Air	0.000	13:45
02 *Subject Test	INT*	13:46
03 Room Air	0.000	13:46

\*Invalid Test

Interferent Detected

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA Cit = INTERFERENT CK Test = OTH

Dr. Lic. = ND/TES989643

Lot No = 14323080A4

Cyl No = 13

Expiration Date = 06/05/2025

Oper No. = 133237County = 08

> Operator Signature ANNA NAREHOOD

Remarks: Interferent check

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004948
Location = TOXL 8164.14.00 09/16
04/01/2024 13:47

Test	AC	Time
01 Diagnostic	OK	13:47
02 Room Air	0.000	13:48
03 *Subject Test	RFI*	13:48
04 Room Air	0.000	13:49

\*Invalid Test Inhibited - RFI

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA Test = OTH Cit = RFI CHECK

Dr. Lic. = ND/TES989643

Lot No = 14323080A4

Cyl No = 13

Expiration Date = 06/05/2025

County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks: RFI Chelk

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = TOXL 04/01/2024

SN 80-004948 8164.14.00 09/16 13:49

### DRY CAL CHECK

Test		AC		Time	
01	Room	Air	0.000	13:50	
02	std.	Gas	0.084	13:50	
03	Room	Air	0.000	13:51	
04	Std.	Gas	0.083	13:51	
05	Room	Air	0.000	13:51	
06	Std.	Gas	0.084	13:52	
07	Room	Air	0.000	13:52	

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Dry cal check- needs calibration adjustment

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/01/2024

SN 80-004948 8164.14.00 09/16 15:34

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:35
02	Std.	Gas	0.079	15:35
03	Room	Air	0.000	15:35
04	Std.	Gas	0.078	15:36
05	Room	Air	0.000	15:36
06	Std.	Gas	0.079	15:37
07	Room	Air	0.000	15:37

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Dry Cal Cheuc- post cal adj

#1-3

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/01/2024

Alcohol Analyzer SN 80-004948 8164.14.00 09/16 15:37

#### DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:38
02	Std.	Gas	0.079	15:38
03	Room	Air	0.000	15:39
04	Std.	Gas	0.079	15:39
05	Room	Air	0.000	15:40
06	Std.	Gas	0.079	15:40
07	Room	Air	0.000	15:40

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature
ANNA NAREHOOD

Remarks:

 $\pi \uparrow \psi$ 

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/01/2024

SN 80-004948 8164.14.00 09/16 15:42

## DRY CAL CHECK

Test		t AC		Time
01	Room	Air	0.000	15:43
02	Std.	Gas	0.079	15:43
03	Room	Air	0.000	15:44
04	Std.	Gas	0.079	15:44
05	Room	Air	0.000	15:45
06	Std.	Gas	0.079	15:45
07	Room	Air	0.000	15:45

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Dry Cal Check #1-9