

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

INTOXILYZER® 8000 ANNUAL INSPECTION

Intoxil	yzer® 8	3000 Serial Number:	80-00 4188	Inspe	ction Location:	TOXL
Α.	Pre-In 1. 2. 3. 4. 5.	Keys Power Cord ☐ Download Data ☐ Upload Operator ☐ Current Location ☐ Battery Check	Yes or No (If Yes, Lo Yes or No Yes or No File		Yes or No	Cyl. #)
		Replaced Sin	nulator O-Ring eath Tube O-Ring		Yes or No Yes or No	
B.	1. 2. 3.	Breath tube heated Date, time and lo Time Zone: (Print test (Level 1 Tank monitor (Le Display: 1 Display and F Gas tank tare If Yes, display	sed and instrument in ed cation code (Level 2, CST or MDT (Time of 1,P). Sign and attached a second	E). Re-set if n test records test record. ator: 175 each other	psi Yes or No G):	zone circled).
C.		Wet Calibration Known Value Sim. Ser #: Lot #: Exp. Date:	ator for the following to Check - Low AC (0.03 AC: 0.02 CM 0.02 CM 0.03 AC:	Level 1,C)	5).	

Intoxilyzer 8000 Annual Inspection

Document ID: 11698 Revision: 2

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Status: Published

Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

Date Approved: 02/29/2024

Page 1 of 2

	2.	Configure simulator for the following test (Level 1,S).					
		Wet Calibration Check - High AC (Level 1,C)					
		Known Value ≥ 0.25 AC: 0.400 AC					
		Sim. Ser #:MP6036					
		Lot #: 202212C					
		Exp. Date: 08Dec 24					
		Results ± 5% AC of known AC					
	3.	Configure dry gas standard for the remaining tests (Level 1,S).					
		Known Value : <u>0.080</u> AC					
		Gas Cylinder Lot #: 14323080A4					
		Cylinder #:					
	4.	Interferent Check (Level 1,B)					
		Known Value: O.IO AC + 0.05% Acetone					
		Sim. Ser #:					
		Lot #: LCS 8					
		Exp. Date:NIA					
		☑ Display reads "Interferent Detect"					
	5.	RFI Check (CMS Mode)					
		Display reads "RFI Detect"					
	6.						
		Test 1 0.081 Test 4 0.081 Test 7 0.081					
		Test 2 0.080 Test 5 0.082 Test 8 0.082					
		Test 3 0.082 Test 6 0.081 Test 9 0.081					
		Average 0.081					
		Results ± 0.005 AC of known AC					
D.	Rema	rks/Maintenance: NIA					
nstru	ment is	acceptable to be used in the field. Yes or No					
f No	state re	eason(s) why:					
,	otato it	sacon(c) miy					
f Yes	chang	ge location code back to A.4. □ NIA					
O2Apra024							
nspector Signature Date							
Strelle Pullschuler 03 Apr 2024							
Reviewer Date							

Intoxilyzer 8000 Annual Inspection

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

Document ID: 11698 Revision: 2

Status: Published

Date Approved: 02/29/2024

Page 2 of 2

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/02/2024 10:31

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 34921080A1 37 Standard Expiration? 02/05/2024 Oper No? 133237

Flow Cal. Date: 04/30/2020 Slope 645 Intercept -637424

IR Calibration Date: 05/06/2020 3um 9um

Oth Coef(*100): -30092 -33449
1st Coef(*100): 293453 143538
2nd Coef(*100): 151 602
H2O adj(mg/l*10k): 425 462

******* Printer Test End *******

Operator Signature ANNA NAREHOOD

Remarks: Print Test

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/02/2024

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 10:37

WET CAL CHECK

Test			AC	Time	
01	Room	Air	0.000	10:38	
02	Std.	Sol.	0.020	10:39	
03	Room	Air	0.000	10:40	
04	Std.	Sol.	0.020	10:41	
05	Room	Air	0.000	10:41	
06	Std.	Sol.	0.021	10:42	
07	Room	Air	0.000	10:42	

08 Sim Temp = 34.0°C

Simul Ser No = MP3002 Std Sol No = 202301E

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: LOW AC Check- 0.020AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/02/2024

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 10:50

WET CAL CHECK

Test			AC		Time
01	Room	Air	0.00	00	10:50
02	Std.	Sol.	0.39	94	10:51
03	Room	Air	0.00	00	10:52
04	Std.	Sol.	0.39	95	10:53
05	Room	Air	0.00	00	10:53
06	Std.	Sol.	0.39	95	10:54
07	Room	Air	0.00	00	10:54

08 Sim Temp = 34.0°C

Simul Ser No = MP6036 Std Sol No = 202212C

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: High AC Check- 0.400AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = TOXL 8164.14.00 09/16
04/02/2024 11:09

Test			AC	Time
01	Room Air		0.000	11:10
02	*Subject	Test	INT*	11:11
03	Room Air		0.000	11:11

*Invalid Test

Interferent Detected

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA

Test = OTH Cit = INTERFERENT CK

Dr. Lic. = ND/TES989643

Lot No = 14323080A4

Cyl No = 13

Expiration Date = 06/05/2025

County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Interferent check

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004188
Location = TOXL	8164.14.00 09/16
04/02/2024	11:12

Test	AC	Time	
01 Diagnostic	OK	11:13	
02 Room Air	0.000	11:14	
03 *Subject Test	RFI*	11:14	
04 Room Air	0.000	11:15	

*Invalid Test Inhibited - RFI

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA Test = OTH Cit = RFI CHECK

Dr. Lic. = ND/TES989643

Lot No = 14323080A4

Cyl No = 13

Expiration Date = 06/05/2025

County = 08

Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks: RFI (Neuk

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/02/2024

SN 80-004188 8164.14.00 09/16 11:15

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	11:16
02	Std.	Gas	0.081	11:16
03	Room	Air	0.000	11:17
04	Std.	Gas	0.080	11:17
05	Room	Air	0.000	11:17
06	Std.	Gas	0.082	11:18
07	Room	Air	0.000	11:18

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

Form 106-I8000

Dry cal cheek- &

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/02/2024

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 11:20

DRY CAL CHECK

Test			AC	Time	
01	Room	Air	0.000	11:21	
02	Std.	Gas	0.081	11:21	
03	Room	Air	0.000	11:21	
04	Std.	Gas	0.082	11:22	
05	Room	Air	0.000	11:22	
06	Std.	Gas	0.081	11:23	
07	Room	Air	0.000	11:23	

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Dry Cal Check #4-6

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
04/02/2024

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 11:23

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	11:24
02	Std.	Gas	0.081	11:24
03	Room	Air	0.000	11:25
04	Std.	Gas	0.082	11:25
05	Room	Air	0.000	11:26
06	Std.	Gas	0.081	11:26
07	Room	Air	0.000	11:27

Lot No = 14323080A4

Cyl No = 13

Exp Date = 06/05/2025

County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Dry Cal Check #7-9