

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 66 91	Instrument Location NEWT	•		
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change				
* Other (Specify) annual Inspection at Crime Laboratory				
Check When Done:	O			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
▼ 6. Rup Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
O. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date 4-3-241		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By		Certified Date 08April2024		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006691
Location = NEWT 8164.14.00 09/16
04/03/2024 19:38

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Yes Display Prelim Rslt? Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? Standard Expiration?

Oper No? 132000

Flow Cal. Date: 08/28/2018 Slope 699 Intercept -660610

IR Calibration Date: 04/03/2024

3um 9um

 0th Coef(*100):
 -12327
 -23187

 1st Coef(*100):
 262018
 135444

 2nd Coef(*100):
 2687
 1323

 H2O adj(mg/l*10k):
 263
 405

****** Printer Test End ********

Operator Signature MARK NEW HOLY

Remarks:

Print Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
04/03/2024

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 19:48

DRY CAL CHECK

Test	AC	Time
01 Room Ad 02 Std. Ga 03 Room Ad 04 Std. Ga 05 Room Ad 06 Std. Ga 07 Room Ad	as 0.079 ir 0.000 as 0.078 ir 0.000 as 0.079	19:49 19:49 19:50 19:50 19:51 19:51
0 / 100M 11.		

Lot No = 14323080A1

Cyl No = 71

Exp Date = 06/05/2025

County = 60

Oper No. = 132000

Operator Signature MARK NEW HOLY

Remarks: ALA Tast

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
04/03/2024

Alcohol Analyzer SN 80-006691 8164.14.00 09/16 19:41

Test		AC	Time
01 Room Air	Test	0.000	19:43
02 *Subject		RFI*	19:46
03 Room Air		0.000	19:46

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST T Sub DOB = 01/12/1999

Sub DOB = 01/12/1999
Sub Sex = Male
Test = DUI
Dr. Lic. = ND/TES991234

Weight = 180 Cit = TEST

Lot No =
Cyl No =

Expiration Date =

County = 60

Oper No. = 132000

Operator Signature MARK NEWHOLY

Remarks:

RFI Test

Form 106-I8000