



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-004935</i>	Instrument Location <i>Linton</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <i>Annual Inspection at Crime Laboratory</i>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RF Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>4/5/24</i>
---	-----------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>[Signature]</i>	Certified Date <b>08April2024</b>
---	--------------------------------------

Intoxilyzer Test Record and Checklist  
NDOPAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004935  
Location = LINT 8164.14.00 09/16  
04/05/2024 17:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	17:25
02 Room Air	RFI*	17:25

\*Invalid Test  
Inhibited - RFI

Lot No = 143230804  
Cyl No = 39  
Exp Date = 06/05/2025  
County = JS  
Oper. No. = 130117

  
Operator Signature  
MARK PEARSON

Remarks:

Form 106-18000

INTOXILYZER 8000  
Instrument Initialization  
16:50 04/05/2024

INTOXILYZER 8000  
Instrument Initialization  
17:07 04/05/2024

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004935  
Location = LINT 8164.14.00 09/16  
04/05/2024 16:57


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004935  
Location = LINT 8164.14.00 09/16  
04/05/2024 17:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:20
02 Std. Gas	0.079	17:21
03 Room Air	0.000	17:21
04 Std. Gas	0.079	17:22
05 Room Air	0.000	17:22
06 Std. Gas	0.080	17:23
07 Room Air	0.000	17:23

Lot No = 14323080A4  
Cyl No = 39  
Exp Date = 06/05/2025  
County = 15 Oper No. = 130117

  
Operator Signature  
MARK PEARSON

Remarks:

Form 106-I8000

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMNPO

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNPOQRSTU

24 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNPOQRSTUUVWX

32 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNPOQRSTUUVWXYZ123456

42 characters per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNPOQRSTUUVWXYZ1234567890ABCDEF

Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt?: Yes  
Display Third Digit?: Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 14323080A4  
Standard Cyl #? 39  
Standard Expiration? 06/05/2025  
Oper No? 130117

Flow Cal. Date: 07/12/2012  
Slope: 655  
Intercept: -510289

IR Calibration Date: 05/20/2015

	3um	9um
0th Coef(*100):	-7871	-19684
1st Coef(*100):	264505	133108
2nd Coef(*100):	2865	1320
H2O Hgmo (*10K):	200	298

\*\*\*\*\* Test End \*\*\*\*\*

  
Operator Signature  
MARK PEARSON

Remarks:

Form 106-I8000