



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION**

**INTOXILYZER® 8000 CALIBRATION ADJUSTMENT**

Intoxilyzer® 8000 Serial Number: 80-00 4947 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No

B. Calibration Adjustment (Level 3,M,C,O)

1.  Autocalibration Printout Attached
  - Max Power Res Value  $\geq 10$
  - Auto Range Res Value  $\geq 4$
2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H <sub>2</sub> O	NA-Milli-Q H <sub>2</sub> O	MP5321
2	0.040	<u>202303H</u>	<u>28Mar25</u>	MP5289
3	0.080	<u>202302B</u>	<u>14Feb25</u>	MP3067
4	0.100	<u>202304A</u>	<u>04Apr25</u>	MP4038
5	0.300	<u>202402C</u>	<u>14Feb26</u>	MP3062

3. 0.080 AC Calibration Gas for H<sub>2</sub>O Adjustment

Lot No. 14323080A4 Cyl No. 13 Exp. Date: 6/5/25

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 929 mbar  
 Adjusted to using barometer 949 mbar  
 Auto Calibration Report printout 949 mbar  
 Barometer Model 10510-922  
 Barometer Serial Number 230307250  
 Barometer Calibration Expiration Date 02May25

5.  Screen displayed "Calibration Success"


6.  Calibration Adjustment Printout Attached

- Solution 1 Avg % Abs  $\leq 0.2500$
- Solution 2-5 REL STD DEV  $\leq 3.000$
- Residual (g/210 L) values for solutions 1 - 5  $\leq 0.0020$  for 3  $\mu$ m and 9  $\mu$ m channels

Dry Gas H<sub>2</sub>O adjustment sum for 3 μm and 9 μm channels within ± 10  
3 μm 3400 (Ave.) + 409 (H<sub>2</sub>O Adj.) = 3809  
9 μm 3371 (Ave.) + 438 (H<sub>2</sub>O Adj.) = 3809

C. Is an Annual Inspection due for this instrument? Yes or No  
If Yes, complete Intoxilyzer 8000 Annual Inspection (Document ID: 11698)  
If No, complete Intoxilyzer 8000 Calibration (Document ID: 11871).

Remarks/Notes: N/A  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

  
\_\_\_\_\_  
Analyst Signature

09 Apr 2024  
\_\_\_\_\_  
Date

  
\_\_\_\_\_  
Reviewer Signature

16 Apr 2024  
\_\_\_\_\_  
Date

TDXL  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-004947  
04/09/2024 14:09:39

Auto Calibration  
Max Power Res Value = 35  
Auto Range Res Value = 16

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004947  
 04/09/2024 14:09:39

Auto Calibration

	<<<<< 3um >>>>>		<<<<< 9um >>>>>	
-----				
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.1280	(0.0070)	0.1920	(0.0250)
Sample #2	0.0950	(0.0950)	0.2040	(0.0150)
Sample #3	0.0690	(0.1580)	0.1850	(0.0550)
Sample #4	0.0930	(0.1770)	0.1800	(0.0620)
Avg % Abs	0.0857	(0.1433)	0.1897	(0.0440)
STD DEV	0.0145	(0.0429)	0.0127	(0.0254)
REL STD DEV	16.889	(29.946)	6.676	(57.631)
-----				
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.7670	(-0.0180)	1.5510	(0.0010)
Sample #2	0.7800	(-0.0070)	1.5400	(0.0400)
Sample #3	0.7810	(0.0160)	1.5500	(0.0460)
Sample #4	0.7850	(0.0200)	1.5490	(0.0420)
Avg % Abs	0.7820	(0.0097)	1.5463	(0.0427)
STD DEV	0.0026	(0.0146)	0.0055	(0.0031)
REL STD DEV	0.338	(150.741)	0.356	(7.160)
-----				
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.4730	(-0.0140)	2.9360	(-0.0020)
Sample #2	1.5140	(-0.0150)	2.9560	(0.0240)
Sample #3	1.5010	(0.0000)	2.9350	(0.0400)
Sample #4	1.5210	(-0.0060)	2.9310	(0.0370)
Avg % Abs	1.5120	(-0.0070)	2.9407	(0.0337)
STD DEV	0.0101	(0.0075)	0.0134	(0.0085)
REL STD DEV	0.671	(107.855)	0.457	(25.262)
-----				
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.8520	(-0.0020)	3.5920	(0.0070)
Sample #2	1.8620	(0.0200)	3.6260	(0.0060)
Sample #3	1.8490	(0.0510)	3.6330	(0.0140)
Sample #4	1.8600	(0.0380)	3.6420	(0.0190)
Avg % Abs	1.8570	(0.0363)	3.6337	(0.0130)
STD DEV	0.0070	(0.0156)	0.0080	(0.0066)
REL STD DEV	0.377	(42.845)	0.221	(50.442)
-----				
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1				
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.2190	(0.0060)	9.8140	(-0.0060)
Sample #2	5.2970	(0.0080)	9.9270	(0.0040)
Sample #3	5.2820	(0.0230)	9.9520	(0.0040)
Sample #4	5.2670	(0.0300)	9.9460	(0.0160)
Avg % Abs	5.2820	(0.0203)	9.9417	(0.0080)
STD DEV	0.0150	(0.0112)	0.0131	(0.0069)
REL STD DEV	0.284	(55.278)	0.131	(86.603)
-----				

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-004947  
 04/09/2024 14:09:39

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
-----			-----		
Zero Order Coef	-210.02			-231.84	
First Order Coef	2642.37			1335.99	
Second Order Coef	19.27			12.47	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	0.000	-0.0003	0.000	0.000	-0.0005
0.040	0.039	0.0008	0.040	0.039	0.0009
0.080	0.080	-0.0004	0.080	0.080	0.0001
0.100	0.100	-0.0000	0.100	0.101	-0.0005
0.300	0.300	0.0000	0.300	0.300	0.0001
-----			-----		

<<<<< 3um >>>>>		<<<<< 9um >>>>>	
-----		-----	
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1			
Sample			
Sample #1	3357.00		3370.00
Sample #2	3417.00		3392.00
Sample #3	3395.00		3376.00
Sample #4	3388.00		3347.00
Avg	3400.0000		3371.6667
STD DEV	15.1327		22.8108
REL STD DEV	0.445		0.677
H2O adjust (mg/l*10k)	409		438

Atmospheric Pressure = 949

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*