



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4947

A. General Setup and Checks:

- 1. Initial start up
  - Power On
  - 5-minute countdown during initial warm-up
  - 1-minute countdown to begin test after standby mode
- 2.  Display Intensity (Ctrl + ↑ or ↓)
- 3. Password Security:
  - Level 1 (Operator)
  - Level 2 (Field Inspector)
  - Level 3 (Lab Personnel)
- 4.  Printer Setup (Level 3,S,P)
- 5.  Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

- 1.  Internal Printer Test (Level 1,P). Sign and attach record.
  - Paper Feed
- 2.  External Printer Test (Level 1,P). Sign and attach record.
- 3. Barcode Reader (Level 1,B) (Tone sounds after each)
  - Operator ID
  - Invalid ID
  - ND License
  - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

- 1.  Continuous Air Blank (Level 2,A)
- 2.  Diagnostic (Level 2,D). Sign and attach test record.
- 3.  Prelim Data Entry (Level 2,E)
- 4. Ser Number/Version Info (Level 2,V):
  - Inst Ser#: 80-004947 DSP Ver: 1362.14
  - Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 13211 s: 0 C: 47.00  
9: 13339 s: 1 B: 45.0

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0→1

2. Function (F) Keys:

F5 – Tone

F9 - Fan

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 927 CF: 1.09

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 214 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 48.0 I: 0.0

B: 45.0 V: 0.0

Cell Temperature (C) between  $47 \pm 0.02$  °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 432 psi Regulator: 475 psi

Display and Regulator  $\pm$  50 psi of each other Yes or No

Gas tank tare necessary? Yes or No

If Yes, display readings after tare (Level 3,M,C,G):

Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

7.  General Setup (Level 3,S,G)

8.  Configure Start Test (Level 3,S,T)

9.  Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1.  Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP5321 Lot #: NIA Exp date: NIA

Known value: 0.000 AC

Results are Zero AC

2.  Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value  $\leq$  0.03 AC: 0.020 AC

Sim ser #: MP3002 Lot #: 202301E Exp date: 17Jan25

Results  $\pm$  0.005 of known AC


3.  Configure wet bath simulator for the following test (Level 1,S)  
 High AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\geq$  0.25 AC: 0.400 AC  
 Sim ser #: MP6036 Lot #: 202212C Exp date: 08Dec24  
 Results  $\pm$  5% of known AC
4.  Configure dry gas standard for the following tests (Level 1,S)  
 Known Value: 0.080 AC  
 Gas Cylinder Lot #: 1432308044 Cylinder #: 13 Exp. Date: 6/5/25
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7344 Lot #: ICS 8 Exp. Date: N/A  
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7111 Lot #: ACE3 Exp. Date: N/A  
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):  
 Sim. Ser #: DR7345 Lot #: MET3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):  
 Sim. Ser #: DR7350 Lot #: IS03 Exp. Date: N/A  
 Display Reads "Interferent Detect"
9.  Ambient Fail (Level 1,B): Display read and test record correct.
10.  Invalid Sample (Level 1,B): Display read and test record correct.
11.  Sequence Aborted (Level 1,B): Display read test record correct.
12.  ABA (Level 1,B): Test record correct.
13.  Calibration check out of tolerance (CMS mode): Display read and test record correct.
14.  RFI Detected (B,C, or CMS mode): Display read and test record correct.
15.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
16.  Deficient Samples (CMS mode): Display read and test record correct.
17.  Improper Sample (CMS mode): Display read and test record correct.
18.  Difference Too Great (CMS mode): Display read and test record correct.
19.  CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):  
 Ave: 0.0821 Std Dev: 0.0007 Rel Std Dev: 0.88%

Maintenance/Repairs/Notes: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Instrument is acceptable to be used in the field.  
If No, state reason(s) why:

Yes or No


\_\_\_\_\_  
\_\_\_\_\_



Inspector Signature

09 Apr 2024

Date



Reviewer

10 Apr 2024

Date



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004947  
Location = TOXL 8164.14.00 09/16  
04/09/2024 09:38

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMN

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMN

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMN


42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMN

Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MemFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 34921080A1  
Standard Cyl #? 40  
Standard Expiration? 02/05/2024  
Oper No? 133237

Flow Cal. Date: 04/07/2020  
Slope 618  
Intercept -565524

IR Calibration Date: 04/07/2020  
3um 9um  
-----  
0th Coef(\*100): -22371 -26151  
1st Coef(\*100): 260822 133081  
2nd Coef(\*100): 1936 1144  
H2O adj(mg/l)\*10K: 495 535

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks: Internal  
Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            09:43

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            34921080A1  
Standard Cyl #?            40  
Standard Expiration?        02/05/2024  
Oper No?                    133237

Flow Cal. Date:            04/07/2020  
Slope                        618  
Intercept                    -565524

IR Calibration Date:        04/07/2020

|                    | 3um    | 9um    |
|--------------------|--------|--------|
| 0th Coef(*100):    | -22371 | -26151 |
| 1st Coef(*100):    | 260822 | 133081 |
| 2nd Coef(*100):    | 1936   | 1144   |
| H2O adj(mg/l*10k): | 495    | 535    |

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks: External Print Test


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            09:46

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 09:47 |
| 02 *Subject Test | REF*  | 09:47 |
| 03 Room Air      | 0.000 | 09:48 |

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = REFUSAL  
Dr. Lic. = ND/DIS821456  
Lot No = 34921080A1  
Cyl No = 40  
Expiration Date = 02/05/2024  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Refusal*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      09:42

DIAGNOSTIC

|                           |      |
|---------------------------|------|
| Voltage/Current Test      | Pass |
| RAM Test                  | Pass |
| EEPROM Checksum Test      | Pass |
| Real Time Clock Test      | Pass |
| DSP Test                  | Pass |
| Analytical Stability Test | Pass |
| Modem Test                | Pass |
| Temp Regulation Test      | Pass |

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      11:05

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 11:06 |
| 02 Std. Sol. | 0.000 | 11:08 |
| 03 Room Air  | 0.000 | 11:08 |
| 04 Std. Sol. | 0.000 | 11:09 |
| 05 Room Air  | 0.000 | 11:09 |
| 06 Std. Sol. | 0.000 | 11:10 |
| 07 Room Air  | 0.000 | 11:11 |


08 Sim Temp = 34.0°C

Simul Ser No = MP5321

Std Sol No = NA

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Zero AC check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      11:27

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 11:28 |
| 02 Std. Sol. | 0.018 | 11:29 |
| 03 Room Air  | 0.000 | 11:29 |
| 04 Std. Sol. | 0.019 | 11:30 |
| 05 Room Air  | 0.000 | 11:30 |
| 06 Std. Sol. | 0.019 | 11:31 |
| 07 Room Air  | 0.000 | 11:32 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3002  
Std Sol No = 202301E  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: *LOW AC CHECK- 0.020AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      11:39

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 11:39 |
| 02 Std. Sol. | 0.392 | 11:40 |
| 03 Room Air  | 0.000 | 11:41 |
| 04 Std. Sol. | 0.397 | 11:42 |
| 05 Room Air  | 0.000 | 11:42 |
| 06 Std. Sol. | 0.397 | 11:43 |
| 07 Room Air  | 0.000 | 11:43 |

08 Sim Temp = 34.0°C

Simul Ser No = MP6036

Std Sol No = 202212C

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: High Ac Check -0.400Ac

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            11:47

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 11:48 |
| 02 *Subject Test | INT*  | 11:49 |
| 03 Room Air      | 0.000 | 11:49 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = ETOH & ACE  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent check.  
0.10% ETOH and 0.05% Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            11:51

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 11:52 |
| 02 *Subject Test | INT*  | 11:53 |
| 03 Room Air      | 0.000 | 11:53 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = ACETONE  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Acetone Interferent check  
0.05%

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      11:54

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 11:55 |
| 02 *Subject Test | INT*  | 11:56 |
| 03 Room Air      | 0.000 | 11:56 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = MEOH  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Methanol* | *Interferent check*  
*0.10%*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      11:58

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 11:59 |
| 02 *Subject Test | INT*  | 12:00 |
| 03 Room Air      | 0.000 | 12:00 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = IPA  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: 0.10% Isopropanol Interferent Check

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            12:04

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | AMB*  | 12:04 |
| 02 Room Air | 0.000 | 12:05 |

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = AMBIENT  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Ambient Fail

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      12:16

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Room Air       | 0.000 | 12:17 |
| 02 Invalid Sample | X.XXX | 12:17 |
| 03 Room Air       | 0.000 | 12:18 |

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = INVALID  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Invalid ~~test~~  
Sample

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      12:19

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 12:19 |
| 02 *Subject Test | ABT*  | 12:20 |
| 03 Room Air      | 0.000 | 12:20 |

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = SEQUENCE ABRTD  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Sequence aborted*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      12:24

| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Room Air       | 0.000 | 12:24 |
| 02 Subject Test 1 | 0.000 | 12:25 |
| 03 Room Air       | 0.000 | 12:25 |
| 04 Reported AC    | 0.000 | 12:25 |

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = ABA  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: ABA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            13:03

| Test              | AC     | Time  |
|-------------------|--------|-------|
| 01 Diagnostic     | OK     | 13:04 |
| 02 Room Air       | 0.000  | 13:04 |
| 03 Subject Test 1 | 0.000  | 13:06 |
| 04 Room Air       | 0.000  | 13:10 |
| 05 Std. Gas       | 0.042* | 13:11 |
| 06 Room Air       | 0.000  | 13:11 |

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = CAL CK OOT  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Calibration check  
out of tolerance.

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            12:33

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 12:33 |
| 02 Reference | RFI*  | 12:34 |
| 03 Room Air  | 0.000 | 12:34 |

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            12:45

| Test             | AC     | Time  |
|------------------|--------|-------|
| 01 Diagnostic    | OK     | 12:46 |
| 02 Room Air      | 0.000  | 12:46 |
| 03 *Subject Test | 0.000* | 12:49 |
| 04 Room Air      | 0.000  | 12:55 |
| 05 Std. Gas      | 0.081  | 12:56 |
| 06 Room Air      | 0.000  | 12:58 |
| 07 *Subject Test | 0.000* | 13:01 |
| 08 Room Air      | 0.000  | 13:02 |


No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                    Cit = DEFICIENT SAMPL  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                    Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            12:07

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | IPS*  | 12:07 |
| 02 Room Air | 0.000 | 12:08 |

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = ~~INVALID~~  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                            Oper No. = 133237

*Improper aen  
4/9/24*



Operator Signature  
ANNA NAREHOOD

Remarks: *Improper sample*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      13:22

| Test              | AC     | Time  |
|-------------------|--------|-------|
| 01 Diagnostic     | OK     | 13:23 |
| 02 Room Air       | 0.000  | 13:23 |
| 03 Subject Test 1 | 0.021* | 13:23 |
| 04 Room Air       | 0.000  | 13:26 |
| 05 Std. Gas       | 0.081  | 13:27 |
| 06 Room Air       | 0.000  | 13:28 |
| 07 Subject Test 2 | 0.000* | 13:30 |
| 08 Room Air       | 0.000  | 13:31 |

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = OTH      Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08      Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004947  
Location = TOXL                    8164.14.00 09/16  
04/09/2024                            13:32


| Test              | AC    | Time  |
|-------------------|-------|-------|
| 01 Diagnostic     | OK    | 13:33 |
| 02 Room Air       | 0.000 | 13:34 |
| 03 Subject Test 1 | 0.000 | 13:34 |
| 04 Room Air       | 0.000 | 13:36 |
| 05 Std. Gas       | 0.081 | 13:37 |
| 06 Room Air       | 0.000 | 13:39 |
| 07 Subject Test 2 | 0.000 | 13:39 |
| 08 Room Air       | 0.000 | 13:40 |
| 09 Reported AC    | 0.000 | 13:34 |

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = CMS  
Dr. Lic. = ND/DIS821456  
Lot No = 14323080A4  
Cyl No = 13  
Expiration Date = 06/05/2025  
County = 08                    Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: CMS TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = TOXL      8164.14.00 09/16  
04/09/2024      13:41

| Test                | AC     | Time  |
|---------------------|--------|-------|
| 01 Room Air         | 0.000  | 13:42 |
| 02 Std. Gas         | 0.080  | 13:42 |
| 03 Room Air         | 0.000  | 13:43 |
| 04 Std. Gas         | 0.081  | 13:43 |
| 05 Room Air         | 0.000  | 13:43 |
| 06 Std. Gas         | 0.082  | 13:44 |
| 07 Room Air         | 0.000  | 13:44 |
| 08 Std. Gas         | 0.082  | 13:45 |
| 09 Room Air         | 0.000  | 13:45 |
| 10 Std. Gas         | 0.082  | 13:45 |
| 11 Room Air         | 0.000  | 13:46 |
| 12 Std. Gas         | 0.082  | 13:46 |
| 13 Room Air         | 0.000  | 13:47 |
| 14 Std. Gas         | 0.083  | 13:47 |
| 15 Room Air         | 0.000  | 13:48 |
| 16 Std. Gas         | 0.082  | 13:48 |
| 17 Room Air         | 0.000  | 13:49 |
| 18 Std. Gas         | 0.082  | 13:49 |
| 19 Room Air         | 0.000  | 13:49 |
| 20 Std. Gas         | 0.082  | 13:50 |
| 21 Room Air         | 0.000  | 13:50 |
| 22 Std. Gas         | 0.082  | 13:51 |
| 23 Room Air         | 0.000  | 13:51 |
| 24 Std. Gas         | 0.083  | 13:52 |
| 25 Room Air         | 0.000  | 13:52 |
| 26 Std. Gas         | 0.082  | 13:52 |
| 27 Room Air         | 0.000  | 13:53 |
| 28 Std. Gas         | 0.082  | 13:53 |
| 29 Room Air         | 0.000  | 13:54 |
| 30 Std. Gas         | 0.082  | 13:54 |
| 31 Room Air         | 0.000  | 13:55 |
| 32 Std. Gas         | 0.083  | 13:55 |
| 33 Room Air         | 0.000  | 13:56 |
| 34 Std. Gas         | 0.083  | 13:56 |
| 35 Room Air         | 0.000  | 13:56 |
| 36 Std. Gas         | 0.082  | 13:57 |
| 37 Room Air         | 0.000  | 13:57 |
| 38 Std. Gas         | 0.083  | 13:58 |
| 39 Room Air         | 0.000  | 13:58 |
| 40 Std. Gas         | 0.082  | 13:58 |
| 41 Room Air         | 0.000  | 13:59 |
| 42 Cal Check Stats: |        |       |
| 43 Average          | 0.0821 |       |
| 44 Std Dev          | 0.0007 |       |
| 45 Rel Std Dev      | 0.88%  |       |



Operator Signature  
ANNA NAREHOOD

Remarks: *Stability Test*

Form 106-I8000