



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 6512

A. General Setup and Checks:

1. Initial start up
  - Power On
  - 5-minute countdown during initial warm-up
  - 1-minute countdown to begin test after standby mode
2.  Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
  - Level 1 (Operator)
  - Level 2 (Field Inspector)
  - Level 3 (Lab Personnel)
4.  Printer Setup (Level 3,S,P)
5.  Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1.  Internal Printer Test (Level 1,P). Sign and attach record.
  - Paper Feed
2.  External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
  - Operator ID
  - Invalid ID
  - ND License
  - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1.  Continuous Air Blank (Level 2,A)
2.  Diagnostic (Level 2,D). Sign and attach test record.
3.  Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):  
Inst Ser#: 80-006512 DSP Ver: 1362.14  
Core Ver: 8167.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12841 s: 0 C: 47.06

9: 13407 s: 1 B: 45

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0→1

2. Function (F) Keys:

F5 – Tone

F9 - Fan

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 950 CF: 1.00

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 221 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 48.5 I: 0.0

B: 44.9 V: 0.0

Cell Temperature (C) between  $47 \pm 0.02$  °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 422 psi Regulator: 450 psi

Display and Regulator  $\pm 50$  psi of each other Yes or No

Gas tank tare necessary? Yes or No

If Yes, display readings after tare (Level 3,M,C,G):

Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

7.  General Setup (Level 3,S,G)

8.  Configure Start Test (Level 3,S,T)

9.  Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1.  Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP5321 Lot #: N/A Exp date: N/A

Known value: 0.000 AC

Results are Zero AC

2.  Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value  $\leq 0.03$  AC: 0.020 AC

Sim ser #: MP3002 Lot #: 202301E Exp date: 17Jan25

Results  $\pm 0.005$  of known AC

3.  Configure wet bath simulator for the following test (Level 1,S)  
 High AC-Wet Calibration Check (Level 1,C):  
 Known Value  $\geq$  0.25 AC: 0.400 AC  
 Sim ser #: MP10036 Lot #: 202212C Exp date: 08DEC 24  
 Results  $\pm$  5% of known AC
4.  Configure dry gas standard for the following tests (Level 1,S)  
 Known Value: 0.080 AC  
 Gas Cylinder Lot #: 28423080A3 Cylinder #: 37 Exp. Date: 11/5/25
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7344 Lot #: ICS8 Exp. Date: N/A  
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):  
 Sim. Ser #: DR7111 Lot #: ACE3 Exp. Date: N/A  
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):  
 Sim. Ser #: DR7345 Lot #: MET3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):  
 Sim. Ser #: DR1350 Lot #: ISO3 Exp. Date: N/A  
 Display Reads "Interferent Detect"
9.  Ambient Fail (Level 1,B): Display read and test record correct.
10.  Invalid Sample (Level 1,B): Display read and test record correct.
11.  Sequence Aborted (Level 1,B): Display read test record correct.
12.  ABA (Level 1,B): Test record correct.
13.  Calibration check out of tolerance (CMS mode): Display read and test record correct.
14.  RFI Detected (B,C, or CMS mode): Display read and test record correct.
15.  Tank Pressure Below Minimum (CMS mode): Display reads correct.
16.  Deficient Samples (CMS mode): Display read and test record correct.
17.  Improper Sample (CMS mode): Display read and test record correct.
18.  Difference Too Great (CMS mode): Display read and test record correct.
19.  CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):  
 Ave: 0.082 Std Dev: 0.0000 Rel Std Dev: 0.0106

Maintenance/Repairs/Notes: NIA

---



---




---


Instrument is acceptable to be used in the field.  
If No, state reason(s) why:

Yes or No

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

  
Inspector Signature

10 Apr 2024  
Date

  
Reviewer

11 Apr. 2024  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-006512  
Location = TOXL 8164.14.00 09/16  
04/10/2024 12:10

\*\*\*\*\* Printer Test \*\*\*\*\*

12 char/line  
abcdefghijkl  
ABCDEFGHIJKL

16 char per line  
abcdefghijklmnop  
ABCDEFGHIJKLMN

21 chars per line  
abcdefghijklmnopqrstu  
ABCDEFGHIJKLMNOPQRSTU

24 characters per line  
abcdefghijklmnopqrstuwx  
ABCDEFGHIJKLMNOPQRSTUWX

32 characters per line  
abcdefghijklmnopqrstuwx123456  
ABCDEFGHIJKLMNOPQRSTUWXYZ123456

42 characters per line  
abcdefghijklmnopqrstuwx1234567890abcdef  
ABCDEFGHIJKLMNOPQRSTUWXYZ1234567890ABCDEF


Current Instrument Setup  
Data Entry Mode: Enabled  
Start Test Sequence: DABACABA  
Display Prelim Rslt? Yes  
Display Third Digit? Yes  
Inhib Printer(Y/N)? No  
Display Volume? No  
Disable On MenFull? Yes  
# of Print Copies? 1  
Select Std (D/W/I)? Dry  
Standard Value? 0.080  
Standard Lot #? 01923080A3  
Standard Cyl #? 25  
Standard Expiration? 02/05/2025  
Oper No? 133237

Flow Cal. Date: 08/20/2015  
Slope 680  
Intercept -660856

IR Calibration Date: 08/21/2015  
3um 9um

0th Coef(*100):	-23051	-19261
1st Coef(*100):	266836	137135
2nd Coef(*100):	3974	1466
H2O adj(mg/l*10k):	785	487

\*\*\*\*\* Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks:

Internal Print  
Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
 North Dakota Model 8000      SN 80-006512  
 Location = TOXL      8164.14.00 09/16  
 04/10/2024      12:11

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890- = |  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890- = |  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:      Enabled  
 Start Test Sequence:      DABACABA  
 Display Prelim Rslt?      Yes  
 Display Third Digit?      Yes  
 Inhib Printer (Y/N)?      No  
 Display Volume?      No  
 Disable On Memfull?      Yes  
 # of Print Copies?      1  
 Select Std (D/W/I)?      Dry  
 Standard Value?      0.080  
 Standard Lot #?      01923080A3  
 Standard Cyl #?      25  
 Standard Expiration?      02/05/2025  
 Oper No?      133237


Flow Cal. Date:      08/20/2015  
 Slope      680  
 Intercept      -660856

IR Calibration Date:      08/21/2015

	3um	9um
-----		

0th Coef (*100):	-23051	-19261
1st Coef (*100):	266836	137135
2nd Coef (*100):	3974	1466
H2O adj (mg/l*10k):	785	487

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
 \_\_\_\_\_  
 Operator Signature  
 ANNA NAREHOOD

Remarks: *External Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                 8164.14.00 09/16  
04/10/2024                        12:11

Test	AC	Time
01 Room Air	0.000	12:12
02 *Subject Test	REF*	12:12
03 Room Air	0.000	12:13

\*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                 Weight = 150  
Test = OTH                        Cit = REFUSAL  
Dr. Lic. = ND/DIS821456  
Lot No = 01923080A3  
Cyl No = 25  
Expiration Date = 02/05/2025  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Refusal*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            12:10

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Diagnostics*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006512  
Location = TOXL      8164.14.00 09/16  
04/10/2024      12:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:31
02 Std. Sol.	0.000	12:31
03 Room Air	0.000	12:32
04 Std. Sol.	0.000	12:33
05 Room Air	0.000	12:33
06 Std. Sol.	0.000	12:34
07 Room Air	0.000	12:34

08 Sim Temp = 34.0°C

Simul Ser No = MP5321

Std Sol No = NA

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Zero Ac Check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006512  
Location = TOXL      8164.14.00 09/16  
04/10/2024      12:36

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:37
02 Std. Sol.	0.020	12:37
03 Room Air	0.000	12:38
04 Std. Sol.	0.020	12:39
05 Room Air	0.000	12:39
06 Std. Sol.	0.021	12:40
07 Room Air	0.000	12:41

08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 202301E

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Low AC check - 0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006512  
Location = TOXL      8164.14.00 09/16  
04/10/2024      12:42

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:43
02 Std. Sol.	0.402	12:44
03 Room Air	0.000	12:45
04 Std. Sol.	0.402	12:45
05 Room Air	0.000	12:46
06 Std. Sol.	0.403	12:47
07 Room Air	0.000	12:47

08 Sim Temp = 34.0°C

Simul Ser No = MP6036  
Std Sol No = 202212C  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: High AC Check - 0.400AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            12:58

Test	AC	Time
01 Room Air	0.000	12:59
02 *Subject Test	INT*	12:59
03 Room Air	0.000	13:00

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = INTER. ACE  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Interferent check*  
*0.05% Acetone*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:00

Test	AC	Time
01 Room Air	0.000	13:01
02 *Subject Test	INT*	13:01
03 Room Air	0.000	13:02

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = INTER. MEOH  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Interferent Check  
0.10% Methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:03

Test	AC	Time
01 Room Air	0.000	13:03
02 *Subject Test	INT*	13:04
03 Room Air	0.000	13:05

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = INTER. ISO  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Interferent check*  
*0.10% Isopropanol*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                 8164.14.00 09/16  
04/10/2024                        13:06

Test	AC	Time
01 Room Air	AMB*	13:07
02 Room Air	0.000	13:08

\*Invalid Test  
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                 Weight = 150  
Test = OTH                        Cit = AMBIENT  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                        Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Ambient Test

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:09

Test	AC	Time
01 Room Air	0.000	13:10
02 Invalid Sample	X.XXX	13:10
03 Room Air	0.000	13:11

\*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = INVALID  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Invalid Sample*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:14

Test	AC	Time
01 Room Air	0.000	13:15
02 *Subject Test	ABT*	13:15
03 Room Air	0.000	13:16

\*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = ABORTED SEQ  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks:

*Sequence Aborted*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:16

Test	AC	Time
01 Room Air	0.000	13:17
02 Subject Test 1	0.000	13:18
03 Room Air	0.000	13:18
04 Reported AC	0.000	13:18

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = ABA  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: ABA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006512  
Location = TOXL                                      8164.14.00 09/16  
04/10/2024    13:19

Test	AC	Time
01 Diagnostic	OK	13:20
02 Room Air	0.000	13:20
03 Subject Test 1	0.000	13:20
04 Room Air	0.000	13:23
05 Std. Gas	0.042*	13:24
06 Room Air	0.000	13:24

\*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                                      Weight = 150  
Test = OTH                                      Cit = CAL CK OOT  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                                      Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Calibration Check  
Out Of Tolerance

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:25

Test	AC	Time
01 Room Air	0.000	13:26
02 *Subject Test	RFI*	13:26
03 Room Air	0.000	13:27

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA NAREHOOD

Remarks: RFI Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                    Alcohol Analyzer  
North Dakota Model 8000                SN 80-006512  
Location = TOXL                            8164.14.00 09/16  
04/10/2024                                    13:27

Test	AC	Time
01 Diagnostic	OK	13:28
02 Room Air	0.000	13:29
03 *Subject Test	0.000*	13:32
04 Room Air	0.000	13:34
05 Std. Gas	0.082	13:36
06 Room Air	0.000	13:37
07 *Subject Test	0.000*	13:40
08 Room Air	0.000	13:43

No RFI Detected

\*Deficient Sample - Value Printed was  
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                            Weight = 150  
Test = OTH                                    Cit = DEFICIENT  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                                    Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the  
instructions displayed by the Intoxilyzer  
in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                 8164.14.00 09/16  
04/10/2024                        13:52

Test	AC	Time
01 Diagnostic	OK	13:53
02 Room Air	IPS*	13:53
03 Room Air	0.000	13:54

\*Improper Sample

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                 Weight = 150  
Test = OTH                        Cit = IMPROPER  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                        Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Improper Sample*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            14:08

Test	AC	Time
01 Diagnostic	OK	14:09
02 Room Air	0.000	14:09
03 Subject Test 1	0.000*	14:11
04 Room Air	0.000	14:13
05 Std. Gas	0.079	14:14
06 Room Air	0.000	14:16
07 Subject Test 2	0.033*	14:16
08 Room Air	0.000	14:17

\*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = DIFF TOO GREAT  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006512  
Location = TOXL                    8164.14.00 09/16  
04/10/2024                            13:43


Test	AC	Time
01 Diagnostic	OK	13:44
02 Room Air	0.000	13:45
03 Subject Test 1	0.000	13:46
04 Room Air	0.000	13:48
05 Std. Gas	0.082	13:49
06 Room Air	0.000	13:50
07 Subject Test 2	0.000	13:51
08 Room Air	0.000	13:52
09 Reported AC	0.000	13:46

Difference OK  
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = OTH                            Cit = IMPROPER CMS *amw*  
Dr. Lic. = ND/DIS821456  
Lot No = 28423080A3  
Cyl No = 37  
Expiration Date = 11/05/2025  
County = 08                            Oper No. = 133237

20 minute waiting period ascertained? Y  
Clean Mouthpiece used and disposed? Y  
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: *CMS Test*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006512  
Location = TOXL                 8164.14.00 09/16  
04/10/2024                        14:22

Test	AC	Time
01 Room Air	0.000	14:22
02 Std. Gas	0.082	14:22
03 Room Air	0.000	14:23
04 Std. Gas	0.082	14:23
05 Room Air	0.000	14:24
06 Std. Gas	0.082	14:24
07 Room Air	0.000	14:25
08 Std. Gas	0.082	14:25
09 Room Air	0.000	14:26
10 Std. Gas	0.082	14:26
11 Room Air	0.000	14:27
12 Std. Gas	0.082	14:27
13 Room Air	0.000	14:28
14 Std. Gas	0.082	14:28
15 Room Air	0.000	14:29
16 Std. Gas	0.082	14:29
17 Room Air	0.000	14:29
18 Std. Gas	0.082	14:30
19 Room Air	0.000	14:30
20 Std. Gas	0.082	14:31
21 Room Air	0.000	14:31
22 Std. Gas	0.082	14:32
23 Room Air	0.000	14:32
24 Std. Gas	0.082	14:33
25 Room Air	0.000	14:33
26 Std. Gas	0.082	14:33
27 Room Air	0.000	14:34
28 Std. Gas	0.082	14:34
29 Room Air	0.000	14:35
30 Std. Gas	0.082	14:35
31 Room Air	0.000	14:36
32 Std. Gas	0.082	14:36
33 Room Air	0.000	14:37
34 Std. Gas	0.082	14:37
35 Room Air	0.000	14:38
36 Std. Gas	0.082	14:38
37 Room Air	0.000	14:39
38 Std. Gas	0.082	14:39
39 Room Air	0.000	14:39
40 Std. Gas	0.082	14:40
41 Room Air	0.000	14:40
42 Cal Check Stats:		
43 Average	0.0820	
44 Std Dev	0.0000	
45 Rel Std Dev	0.01%	

Operator Signature  
ANNA NAREHOOD

Remarks: *Stability Test*

Form 106-I8000