



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-004206</b>	Instrument Location <b>Cass County Jail</b>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <b>Annual Inspection at Crime Laboratory</b>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): None
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <b>Kathryn M. Allen #3671</b>	Date <b>4-11-24</b>
--	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <b>Ane Nahl</b>	Certified Date <b>16April2024</b>
--	--------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = CASS                    8164.14.00 09/16  
04/11/2024                            18:19

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            28423080A3  
Standard Cyl #?            10  
Standard Expiration?       11/05/2025  
Oper No?                    090801

Flow Cal. Date:            06/10/2009  
Slope                        690  
Intercept                  -645883

IR Calibration Date:        06/21/2023  
                                  3um                            9um

-----  
0th Coef(\*100):            -14868                        -20907  
1st Coef(\*100):            263667                        130872  
2nd Coef(\*100):            1957                            1123  
H2O adj(mg/l\*10k):        471                            620

\*\*\*\*\* Printer Test End \*\*\*\*\*

#3071  
Operator Signature  
KATHRYN ALLEN

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = CASS      8164.14.00 09/16  
04/11/2024      18:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:20
02 Std. Gas	0.080	18:20
03 Room Air	0.000	18:20
04 Std. Gas	0.081	18:21
05 Room Air	0.000	18:21
06 Std. Gas	0.081	18:22
07 Room Air	0.000	18:22

Lot No = 28423080A3  
Cyl No = 10  
Exp Date = 11/05/2025  
County = 09      Oper No. = 090801

  
Operator Signature  
KATHRYN ALLEN

Remarks:

All Good!!  
Form 106-18000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = CASS      8164.14.00 09/16  
04/11/2024      18:24

Test	AC	Time
01 Room Air	0.000	18:25
02 *Subject Test	RFI*	18:25
03 Room Air	RFI*	18:25

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST TEST  
Sub DOB = 01/01/1901  
Sub Sex = Female      Weight = 100  
Test = DUI      Cit = 01000000000000  
Dr. Lic. = ND/NA  
Lot No = 28423080A3  
Cyl No = 10  
Expiration Date = 11/05/2025  
County = 09      Oper No. = 090801

  
Operator Signature  
KATHRYN ALLEN

Remarks:

RFI test

Form 106-I8000

OK!