

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006.502	Instrument Location	Down stair
Reason for Install/Repair	•	
│	fter Location Change	
X Other (Specify) ONNIA! INSPECTION OF	t Chime Labora	tory
CALIFORNIA TO THE PART OF THE		
Check Witen Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Loc	ation; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S	i).	
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode	or Level 1, Function B or C	Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SF	N50496, Form 120-G) and F	Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 12	0-G) at the Intoxilyzer® Loca	ation at the Agency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	Intoxilyzer® 8000 Installatio	n and Repair Checkout (SFN59281, Form
Field Inspector Signature		Date / / 211
A Comment	168	4-76-09
and the same line Only		
Crime Laboratory Use Only  This installation has been reviewed and the instrument is approved	to be used for the analysis o	f breath to determine alcohol concentration from
This installation has been reviewed and the instrument is approved the date the Field Inspector performed the installation. This record county of Burleigh, North Dakota, is certified to be a true and correct		by Ochloral, Othino Laboratory
Reviewed/Certified By	10	Certified Date 17April2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006502
Location = STUT	8164.14.00 09/16
04/16/2024	19:20

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= | ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Ves

Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	14323080A1
Standard Cyl #?	20
Standard Expiration?	06/05/2025

Flow Cal. Date:	08/19/2015
Slope	674
Intercept	-654491

IR Calibration Date	e:	08/19/2015
	3um	9um
_ <del> </del>	<del>-</del>	·
0th Coef(*100):	-24446	-18103
1st Coef(*100):	267787	135417
2nd Coef(*100):	3387	1610
H2O adj(ng/l*10k):	740	513

132539

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
JOSHUA FASTNACHT

Remarks:

Oper No?

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = STUT
04/16/2024

Alcohol Analyzer SN 80-006502 8164.14.00 09/16 19:27

## DRY CAL CHECK

_T	est		AC	Time
01	Room	Air	0.000	19:28
02	Std.	Gas	0.082	19:29
03	Room	Air	0.000	19:29
04	Std.	Gas	0.081	19:30
05	Room	Air	0.000	19:30
06	Std.	Gas	0.081	19:31
07	Room	Air	0.000	19:31

Lot No = 14323080A1

Cyl No = 20

Exp Date = 06/05/2025

County = 47

Oper No. = 132539

Operator Signature JOSHUA FASTNACHT

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = STUT
04/16/2024

Alcohol Analyzer SN 80-006502 8164.14.00 09/16 19:31

## DRY CAL CHECK

Test	AC.	Time
01 Room Air	RFI*	19:32
02 Room Air	0.000	19:33

\*Invalid Test Inhibited - RFI

Lot No = 14323080A1

Cyl No = 20

Exp Date = 06/05/2025

County = 47

Oper No. = 132539

Operator Signature JOSHUA FASTNACHT

Remarks:

Form 106-I8000