

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006684	Instrument Location West Forgo Police Dot. (WFAD)			
Reason for Install/Repair	Tango Island			
☐ Install After Receiving From Crime Laboratory ☐ Install Aft	er Location Change			
(Specify) Annual Inspection at	Crime Laboratory			
Check When Done:	O			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
💢 3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
☑ 6. Run Tests:				
☑ B. ACA Test (Level 1, Function C).				
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form			
Field Inspector Signature	Date 4/17/2024			
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on a County of Burleigh, North Dakota, is certified to be a true and correct co				
Reviewed/Certified By Au Au Au Au Au Au Au Au Au A	Certified Date 18April2024			

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006684 Location = WFPD 8164.14.00 09/16 04/17/2024

******* Printer Test ********

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA

Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080

Standard Lot #? Standard Cyl #? 14323080A1

25

Standard Expiration? 06/05/2025 Oper No? 121019

Flow Cal. Date: 09/06/2016 Slope 709 Intercept -725441

IR Calibration Date: 09/06/2016

3um Oth Coef(*100): -16013 1st Coef(*100): 256431 2nd Coef(*100): 2904 -23467 135731 1268 H20 adj (mg/l*10k): 553514

****** Printer Test End ********

Operator Signature RYAN HASKELL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WFPD
04/17/2024

Alcohol Analyzer SN 80-006684 8164.14.00 09/16 21:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:30
02 Std. Gas	0.082	21:30
03 Room Air	0.000	21:31
04 Std. Gas	0.082	21:31
05 Room Air	0.000	21:32
06 Std. Gas	0.081	21:32
07 Room Air	0.000	21:33

Lot No = 14323080A1

Cyl No = 25

Exp Date = 06/05/2025

County = 09

Oper No. = 121019

Operator Signature RYAN HASKELL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

SN 80-006684

21:34

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = WFPD 8164.14.00 09/16 04/17/2024

Test			AC	Time
01	Room	Air	RFI*	21:36
02	Room	Air	RFI*	21:36

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, INSTALL NA

Sub DOB = 09/13/1982

Sub Sex = MaleWeight = 170 Test = OTH Cit = NA

Dr. Lic. = ND/NALot No = 14323080A1

Cyl No = 25

Expiration Date = 06/05/2025

County = 09Oper No. = 121019

> Operator Signature RYAN HASKELL

Remarks:

Form 106-I8000