



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 2667

A. General Setup and Checks:

1. Initial start up
 - Power On
 - 5-minute countdown during initial warm-up
 - 1-minute countdown to begin test after standby mode
2. Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - Level 1 (Operator)
 - Level 2 (Field Inspector)
 - Level 3 (Lab Personnel)
4. Printer Setup (Level 3,S,P)
5. Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1. Internal Printer Test (Level 1,P). Sign and attach record.
 - Paper Feed
2. External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
 - Operator ID
 - Invalid ID
 - ND License
 - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1. Continuous Air Blank (Level 2,A)
2. Diagnostic (Level 2,D). Sign and attach test record.
3. Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-002667 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12521 s: 0 C: 47.00
9: 12496 s: 1 B: 45

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0>1

2. Function (F) Keys:

F5 – Tone

F9 - Fan (does not turn on. classroom only inst.)

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 955 CF: 1.05

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 71 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 47.5 I: 00

B: 44.9 V: 0.0

Cell Temperature (C) between 47 ± 0.02 °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 518 psi Regulator: 515 psi

Display and Regulator ± 50 psi of each other Yes or No

Gas tank tare necessary? Yes or No

If Yes, display readings after tare (Level 3,M,C,G):

Display: _____ psi Regulator: _____ psi

7. General Setup (Level 3,S,G)

8. Configure Start Test (Level 3,S,T)

9. Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1. Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP5321 Lot #: NA Exp date: N/A

Known value: 0.000 AC

Results are Zero AC

2. Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value ≤ 0.03 AC: 0.020 AC

Sim ser #: MP3002 Lot #: 202301E Exp date: 17Jan25

Results ± 0.005 of known AC

3. Configure wet bath simulator for the following test (Level 1,S)
 High AC-Wet Calibration Check (Level 1,C):
 Known Value ≥ 0.25 AC: 0.400 AC
 Sim ser #: MP6036 Lot #: 202212C Exp date: 11/08Dec24
 Results $\pm 5\%$ of known AC
4. Configure dry gas standard for the following tests (Level 1,S)
 Known Value: 0.080 AC
 Gas Cylinder Lot #: 28423080A³ Cylinder #: 11 Exp. Date: 11/5/25
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7344 Lot #: 1CS 8 Exp. Date: NIA
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7111 Lot #: ACE 3 Exp. Date: NIA
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):
 Sim. Ser #: DR7345 Lot #: MET3 Exp. Date: NIA
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):
 Sim. Ser #: DR7350 Lot #: ISO 3 Exp. Date: NIA
 Display Reads "Interferent Detect"
9. Ambient Fail (Level 1,B): Display read and test record correct.
10. Invalid Sample (Level 1,B): Display read and test record correct.
11. Sequence Aborted (Level 1,B): Display read test record correct.
12. ABA (Level 1,B): Test record correct.
13. Calibration check out of tolerance (CMS mode): Display read and test record correct.
14. RFI Detected (B,C, or CMS mode): Display read and test record correct.
15. Tank Pressure Below Minimum (CMS mode): Display reads correct.
16. Deficient Samples (CMS mode): Display read and test record correct.
17. Improper Sample (CMS mode): Display read and test record correct.
18. Difference Too Great (CMS mode): Display read and test record correct.
19. CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):
 Ave: 0.0815 Std Dev: 0.0008 Rel Std Dev: 0.94%


Maintenance/Repairs/Notes: _____

Instrument is acceptable to be used in the field.

Yes or No

If No, state reason(s) why:


per prior annual inspections, this instrument is
classroom only due to the fan motor not turning
on.



Inspector Signature

18 Apr 2024

Date



Reviewer

18 Apr 2024

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 13:59

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMN

24 characters per line
abcdefghijklmnopqrstuwx
ABCDEFGHIJKLMN

32 characters per line
abcdefghijklmnopqrstuwx123456
ABCDEFGHIJKLMN

42 characters per line
abcdefghijklmnopqrstuwx1234567890abcdef
ABCDEFGHIJKLMN


Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.020
Standard Solution #? 202301E
Sim Serial #? MP3002
Oper No? 133237

Flow Cal. Date: 04/19/2021
Slope 615
Intercept 113911

IR Calibration Date: 05/08/2017
3um 9um

0th Coef(*100): -16077 -8130
1st Coef(*100): 268518 134907
2nd Coef(*100): 2115 1175
H2O adj(mg/l*10K): 380 220

***** Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: Internal Print
Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:02

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Wet
Standard Value? 0.020
Standard Solution #? 202301E
Sim Serial #? MP3002
Oper No? 133237

Flow Cal. Date: 04/19/2021
Slope 615
Intercept 113911

IR Calibration Date: 05/08/2017
 3um 9um

0th Coef(*100): -16077 -8130
1st Coef(*100): 268518 134907
2nd Coef(*100): 2115 1175
H2O adj(mg/l*10k): 380 220

***** Printer Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: External Print Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:08

Test	AC	Time
01 Room Air	0.000	14:09
02 *Subject Test	REF*	14:09
03 Room Air	0.000	14:10

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = REFUSAL
Dr. Lic. = ND/DIS821456
Simul Ser No = MP5321
Std Sol No = NA Oper No. = 133237
County = 08



Operator Signature
ANNA NAREHOOD

Remarks: Refusal

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:00

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:03

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:04
02 Std. Sol.	0.000	14:05
03 Room Air	0.000	14:05
04 Std. Sol.	0.000	14:06
05 Room Air	0.000	14:06
06 Std. Sol.	0.000	14:07
07 Room Air	0.000	14:07

08 Sim Temp = 34.0°C

Simul Ser No = MP5321
Std Sol No = NA Oper No. = 133237
County = 08



Operator Signature
ANNA NAREHOOD

Remarks: *zero Ac check*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:13

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:13
02 Std. Sol.	0.022	14:14
03 Room Air	0.000	14:15
04 Std. Sol.	0.022	14:16
05 Room Air	0.000	14:16
06 Std. Sol.	0.022	14:17
07 Room Air	0.000	14:17


08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 202301E

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *LOW AC Check - 0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:28

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:29
02 Std. Sol.	0.397	14:29
03 Room Air	0.000	14:30
04 Std. Sol.	0.399	14:31
05 Room Air	0.000	14:31
06 Std. Sol.	0.400	14:32
07 Room Air	0.000	14:32


08 Sim Temp = 34.0°C

Simul Ser No = MP6036

Std Sol No = 202212C

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: High Ac check - 0.400Ac

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:35

Test	AC	Time
01 Room Air	0.000	14:36
02 *Subject Test	INT*	14:37
03 Room Air	0.000	14:37

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ICS
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Ethanol + 0.05% Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:42

Test	AC	Time
01 Room Air	0.000	14:43
02 *Subject Test	INT*	14:43
03 Room Air	0.000	14:44

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ACE
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Interferent check*
0.05% Acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:39

Test	AC	Time
01 Room Air	0.000	14:40
02 *Subject Test	INT*	14:41
03 Room Air	0.000	14:42

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982 Weight = 150
Sub Sex = Female Cit = INTERFERENT MET
Test = OTH
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Methanol

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:37

Test	AC	Time
01 Room Air	0.000	14:38
02 *Subject Test	INT*	14:39
03 Room Air	0.000	14:39

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ISO
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Interferent check*
0.10% Isopropanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:45

Test	AC	Time
01 Room Air	AMB*	14:46
02 Room Air	0.000	14:46

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = AMBIENT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Ambient Fail

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:52

Test	AC	Time
01 Room Air	0.000	14:52
02 Invalid Sample	X.XXX	14:53
03 Room Air	0.000	14:53

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INVALID
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Invalid Sample

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:54

Test	AC	Time
01 Room Air	0.000	14:55
02 *Subject Test	ABT*	14:55
03 Room Air	0.000	14:55

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ABORTED SEQ
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Sequence aborted*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 15:27

Test	AC	Time
01 Diagnostic	OK	15:28
02 Room Air	0.000	15:28
03 Subject Test 1	0.000	15:29
04 Room Air	0.000	15:31
05 Std. Gas	0.000*	15:32
06 Room Air	0.000	15:33

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CAL CHECK OOT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: *cal check out of tolerance*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:55

Test	AC	Time
01 Room Air	0.000	14:56
02 *Subject Test	RFI*	14:57
03 Room Air	0.000	14:58

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ~~ABA~~ RFI AEN
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 14:58

Test	AC	Time
01 Diagnostic	OK	14:59
02 Room Air	0.000	14:59
03 *Subject Test	0.000*	15:03
04 Room Air	0.000	15:05
05 Std. Gas	0.081	15:07
06 Room Air	0.000	15:08
07 *Subject Test	0.000*	15:11
08 Room Air	0.000	15:12

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ~~DIFF TOO GREAT~~
Dr. Lic. = ND/DIS821456 *deficient AEN*
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.

Operator Signature
ANNA NAREHOOD

Remarks: *Deficient Samples*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 15:25

Test	AC	Time
01 Diagnostic	OK	15:26
02 Room Air	IPS*	15:26
03 Room Air	0.000	15:27

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = IMPROPER
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Improper Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 15:33

Test	AC	Time
01 Diagnostic	OK	15:33
02 Room Air	0.000	15:34
03 Subject Test 1	0.043*	15:34
04 Room Air	0.000	15:37
05 Std. Gas	0.081	15:38
06 Room Air	0.000	15:39
07 Subject Test 2	0.000*	15:39
08 Room Air	0.000	15:40

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 11
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: *Difference Too Great*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002667
Location = TOXL 8164.14.00 09/16
04/18/2024 15:41

Test	AC	Time
01 Room Air	0.000	15:41
02 Std. Gas	0.081	15:41
03 Room Air	0.000	15:42
04 Std. Gas	0.081	15:42
05 Room Air	0.000	15:43
06 Std. Gas	0.081	15:43
07 Room Air	0.000	15:44
08 Std. Gas	0.083	15:44
09 Room Air	0.000	15:44
10 Std. Gas	0.081	15:45
11 Room Air	0.000	15:45
12 Std. Gas	0.081	15:46
13 Room Air	0.000	15:46
14 Std. Gas	0.081	15:47
15 Room Air	0.000	15:47
16 Std. Gas	0.081	15:47
17 Room Air	0.000	15:48
18 Std. Gas	0.081	15:48
19 Room Air	0.000	15:49
20 Std. Gas	0.081	15:49
21 Room Air	0.000	15:50
22 Std. Gas	0.081	15:50
23 Room Air	0.000	15:50
24 Std. Gas	0.081	15:51
25 Room Air	0.000	15:51
26 Std. Gas	0.083	15:52
27 Room Air	0.000	15:52
28 Std. Gas	0.082	15:52
29 Room Air	0.000	15:53
30 Std. Gas	0.082	15:53
31 Room Air	0.000	15:54
32 Std. Gas	0.082	15:54
33 Room Air	0.000	15:55
34 Std. Gas	0.083	15:55
35 Room Air	0.000	15:55
36 Std. Gas	0.081	15:56
37 Room Air	0.000	15:56
38 Std. Gas	0.081	15:57
39 Room Air	0.000	15:57
40 Std. Gas	0.082	15:57
41 Room Air	0.000	15:58
42 Cal Check Stats:		
43 Average	0.0815	
44 Std Dev	0.0008	
45 Rel Std Dev	0.94%	



Operator Signature
ANNA NAREHOOD

Remarks:

Stability Test

Form 106-I8000