

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

80-064942	strument Location Cass County Jail Ma	ichine#2		
Reason for Install/Repair	J .			
☐ Install After Receiving From Crime Laboratory ☐ Install After L	ocation Change			
(Specify) Ohnual Inspection at C	rime Laboratory			
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.				
☑ 3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):	int. Needed			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5049)		oxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxil 104-G), Print Test, ACA Test, and RFI Test.	lyzer® 8000 Installation and Repair Ch	eckout (SFN59281, Form		
Field Inspector Signature  Palfugn M. allu #367		Date 4-22-34		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be used the date the Field Inspector performed the installation. This record on file a County of Burleigh, North Dakota, is certified to be a true and correct copy	at the Office of Attorney General, Crime			
Reviewed/Certified By		Certified Date 29April2024		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = CASS 8164.14.00 09/16
04/22/2024 17:40

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? # of Print Copies? Yes 1 Dry 0.080 Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? 01923080A3 Standard Cyl #? 23 Standard Expiration? 02/05/2025 090801 Oper No?

Flow Cal. Date: 05/25/2018 Slope 626 Intercept -548004

IR Calibration Date: 04/09/2024

3um 9um

Oth Coef(\*100): -22609 -23575
1st Coef(\*100): 260445 132718
2nd Coef(\*100): 2885 1338
H2O adj(mg/l\*10k): 492 488

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

KATHRYN ALLEN

Remarks:

2000 test!

erator Signature

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = CASS 8164.14.00 09/16 04/22/2024

18:06

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	18:07
02	Std.	Gas	0.080	18:07
03	Room	Air	0.000	18:08
04	Std.	Gas	0.080	18:08
05	Room	Air	0.000	18:09
06	Std.	Gas	0.080	18:09
07	Room	Air	0.000	18:10

Lot No = 01923080A3

Cyl No = 23

Exp Date = 02/05/2025

County = 09

Oper No. = 090801

Remarks:

Form \$6-18000

tor Signature

Good test All in Normal

ATHRYN ALLEN

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = CASS 8164.14.00 09/16
04/22/2024 18:10

Test	AC	Time
01 Room Air 02 *Subject Test	0.000 (RFI*	18:11 18:12
03 Room Air	0.000	18:12

\*Invalid Test Inhibited - RFI

Sub Name = TEST AFTER INSP, TEST ATER INSP TEST

Sub DOB = 01/01/1901

Sub Sex = Female

Weight = 225

Test = DUI

Cit = 010100000000000

Dr. Lic. = ND/NA

Lot No = 01923080A3

Cyl No = 23

Expiration Date = 02/05/2025

County = 09

Oper No. = 090801

Operator Signature KATHRYN ALLEN

Remarks:

Form 106-I8000